



User Manual

Pulse Oximeter

MD1913 FSD20



FLEMING MEDICAL LTD
21/F, Chinachem Tower,
34-37 Connaught Road Central,
HONG KONG

FLEMING MEDICAL LTD
Limerick, V94 K6W4, Ireland
1800 307777 (IRE) 0800 7314315 (UK)
WWW.FLEMINGMEDICAL.IE

Rev.:07072020
File Ver:2.1
Release Date: June 2020

Instructions to User

Dear users, thank you very much for purchasing the Pulse Oximeter.
This manual contains the instructions necessary to operate the product safely and in accordance with its function and intended use.
Read this manual carefully before using the fingertip pulse oximeter.
This product is reusable medical device. The using life is 2 years.
The device is prescription use device.

1 Safety

- 1.1.1 Contraindications
Do NOT use oximeter in a magnetic resonance (MR or CT) environment.
1.1.2 Warnings
Keep the oximeter away from young children. Small parts such as the battery door, battery and lanyard may trigger choking hazards.
1.1.3 Cautions
Do not use oximeter in the present of flammable anesthetics.
The oximeter needs to be used according to information provided in the user manual.
The equipment is NOT intended for neonates and infants.
Do not use a damaged oximeter which may affect measurement performance.
Do not place the oximeter on the same hand/arm when using a blood pressure cuff or monitor.
Do not use the oximeter for more than 5 minutes without relocating the device to another finger.
Do not place the oximeter on edema or fragile tissues.
Do not use the oximeter as the only basis for making medical decisions, it is intended only to be used as additional information that you can give to your licensed health care professional.
Do not use the oximeter in high frequency environment such as around electrosurgical equipment.
Do not place the oximeter in liquid.
Follow local disposal and recycling laws for the oximeter and its components, including the battery.
Do not stare at the light (the infrared is invisible) which is emitted from the oximeter and can be harmful to the eyes.
For clinical limitations and contraindications, please carefully review the medical literature.
The equipment is just a clinical diagnosis of auxiliary equipment. The physiological data displayed on the equipment are for reference only and can not be directly used for diagnostic interpretation.
Not suitable for the users with arrhythmia / heart failure / Low perfusion (PI<0.3) / finger trembling.
Not suitable for the users with large finger size or exceeding pulse oximeter's finger measurement cavity size.
Please do not use the thumb or tail finger to measure.
Discomfort or pain may appear if using equipment ceaselessly, especially for microcirculation barrier patients.
It is recommended that the equipment should not be used on the same finger for more than 5 minutes.
The oximeter is designed to measure the percentage of arterial oxygen saturation of functional hemoglobin. Any of the following conditions may reduce the performance of the oximeter.
Flicking or very bright light
Moisture in the oximeter
The individual weight less than 20kg
Weak pulse quality(low perfusion)
Venous pulsations
Low hemoglobin
Cardiogreen and other intravascular dyes
Carboxyhemoglobin
Methemoglobin
Dysfunctional hemoglobin
Artificial nails or fingernail polish
The finger is too cold
Patients with abnormal circulation of finger endings caused by copd

2 The Basics

The pulse oxygen saturation is the percentage of HbO2 in the total Hb in the blood, so-called the O2 concentration in the blood. It is an important bio-parameter for respiration. A number of diseases relating to respiratory system may cause the decrease of SpO2 in the blood, furthermore, some other causes such as the malfunction of human body's self-adjustment, damages during surgery, and the injuries caused by some medical checkup would also lead to the difficulty of oxygen supply in human body. And the corresponding symptoms would appear as a consequence, such as vertigo, impotence, vomit etc. serious symptoms might bring danger to human's life. Therefore, prompt information of patients SpO2 is of great help for the doctor to discover the potential danger, and is of great importance in the clinical medical field.

2.1 Principle
Principle of the oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristic of Reductive hemoglobin(Hb)and oxyhemoglobin(HbO2)in glow & near-infrared zones. Operation Principle of the device is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning Recording Technology, So that two beams of different wavelength of lights can be focused onto human nail tip through perspective clamp finger-type sensor. Then a measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.




Figure 1. Oximeter schematic diagram

2.2 Introduction
2.2.1 Intended Use
The Pulse Oximeter is portable, convenient, non-invasive device, used to monitor arterial hemoglobin oxygen saturation (SPO2) and pulse rate. The personal application are adult patients (weight: >30kg) and pediatric patients (weight: 20-30kg). We recommend index finger, middle finger and ring finger in suitable position for the monitor. It's intended for spot-checking or attended-care monitoring in Home Health Care and Medical Facility.

WARNING:
This pulse oximeter is intended for use only by clinical professionals or under their guidance. It must only be used by persons who have received adequate training in its use. Anyone unauthorized or untrained must not perform any operation on it.

CAUTIONS:
This pulse oximeter is intended for use in hospital,clinical institution,healthcare community.
The pulse oximeter is NOT designed for newborns and infants. For adults and children, it recommended that the finger thickness should between 8-25.4mm.

NOTES
The probe is the hole in the middle of the equipment to which the finger is positioned over.
The probe is the Applied Part of the equipment.

2.2.2 Features
The pulse oximeter is small in volume, light in weight and easy to carry.
One button and easy to operate.
There are three modes: power off, sleep and measure.
Automatically turning into sleep mode within 8 seconds after there is no signal.

NOTES
Press the operating button to activate oximeter (measure mode) from sleep mode.
Power off after removing the batteries.

2.3 Front View




Figure 2. Front View of MD1913

2.4 Functions

Table with 2 columns: Function, MD1913. Rows include Display (OLED), Spo2 parameter measurement (Y), Pulse rate parameter measurement (Y), Bar graph display (Y), Battery display (Y), Automatically enters sleep mode (Y), Pulse waveform display (Y), and Four direction display (Y).

2.5 Symbols

Table with 8 columns: Symbol, Definition, Symbol, Definition, Symbol, Definition, Symbol, Definition. Includes symbols for % SpO2, PR, IPX2, Serial Number, Battery positive electrode, Battery cathode electrode, Date of Manufacturer, and CE mark.

3 Battery Installation

- A. Put the two AAA batteries into battery compartment in correct polarities.
- B. Push the battery cover horizontally along the arrow shown as right.

WARNINGS:

- Battery polarities should be correctly installed, otherwise, damage may be caused to the equipment.
- Please remove the batteries if the oximeter will be stored for more than 30 days.
- Please remove the batteries if you want turn off the oximeter otherwise it is always in power state.
- Battery may leak or explode if used or disposed off improperly




Figure 3. MD1913 Battery Installation

4 Operating Guide

- 4.1 Application Method
A. Remove the battery cover, and insert the two AAA batteries following polarity markings indicated inside of the battery compartment, then reposition the cover;
B. Hold the oximeter with the display facing toward you, slide your finger into the opening probe of the device, as shown below (Figure 4), until the fingertip touches the built-in stop guide. For best results, make sure the finger is centered within the finger guide.
C. Press the button to activate the oximeter from sleep mode, and then measurement interface will appear in 3 seconds.
D. The oximeter will turn into sleep mode automatically within 8 seconds after the finger is removed from the probe.




Figure 4. Operation Guide

- 4.2 Attention for Operation
A. Ring finger, middle finger and index finger are recommended as suitable monitor position.
B. Excessive or rapid movement may affect measurement.
C. Improper sensor placement may affect the measurement accuracy.
D. The oximeter can be reused after cleaning and disinfection.
E. The measurement works best when the oximeter and the heart are at the same level.
F. The plethysmogram can be used as a pulse intensity indicator. The displayed parameters might be unreliable with the disorderly plethysmogram.
G. The displayed parameters will show invalid indicator as '---' if pulse intensity is very low.
H. The displayed parameters will show invalid indicator as '---' if oximeter fault occur.
I. The maximum continuous test time does not exceed 5 minutes.

5 Specifications

- 5.1 Classification
Type of protection against electric shock: II (Internally powered equipment)
Degree of protection against electric shock: Type BF-Applied part (non-defibrillation proof)
Operating mode: Spot checking
Degree of protection against hazards of explosion: Ordinary equipment: Not protected
Equipment type: Fingertip oximeter

- 5.2 Measurement Specifications
Spo2 declared accuracy
Range(σ\*) 70%~100%: ±2%
0% ~ 69%: unspecified
Resolution 1%
Update Period 1s
Averaging Time 8s
PR declared accuracy
Range(σ\*) 25~250 : ±3 BPM (Beats per Minute)
Resolution 1 bpm
Update Period 1s
Averaging Time 8s

- 5.3 Power Requirements
Specification of batteries Two 1.5V(AAA)
Operating voltage DC 2.5V - 3V

- 5.4 Environmental Specifications
Temperature
Operating +41° to +104°F / 5° to +40°C
Storage/Transportation -4° to +140°F / -20° to +60°C
Humidity
Operating 10~95%, noncondensing
Storage/Transportation 10~95%, noncondensing
Atmosphere Pressure
Operating 70~106kpa
Storage/Transportation 50~107.4kpa

- 5.5 Physical Specifications
Width\*Height\*Depth About 33x36x58mm
Weight About 60 g (including the batteries)

5.6 Display

Table with 2 columns: Display type, Display content. Rows include OLED, 0.96", 128 x 64 pixel and SpO2%, Pulse rate, Battery indicator, Bar graph, Pulse waveform, Heart beat mark.

5.7 LED Wavelengths

Table with 3 columns: Wavelength, Radiant Power. Rows include RED 660±6n m 1.8m W and IR 905±10n m 2.0m W.


6 Technical Description

The below table shows statistic conclusion of an invasive controlled desaturation study which guided by 'ISO80601-2-61, Annex EE, Guideline for evaluating and documenting SpO2 Accuracy in human subjects'. The statistic result displayed the accuracy distribution between the range of 70% ~ 100%, which may helpful to user.

Bias Analysis SpO2-MD1913 Pulse Oximeter	SaO2-Radiometer ABL800 FLEX-CO-Oximeter			
	70 – 80 (%)	80 – 90 (%)	90 – 100 (%)	70 – 100 (%)
Mean Bias (B <sub>i</sub> )	1.94	1.45	0.89	1.4
Precision (S <sub>res</sub> )	2	1.55	0.98	1.53
Accuracy (A <sub>rm</sub> )	1.98	1.53	0.96	1.52

The below is the Bland-Altman graphical plot of samples from invasive controlled desaturation study.

Bland-Altman Graph for SaO2-SpO2



7 Maintenance,Cleaning,Disinfection

- 7.1 Maintenance**  
The equipment's design life expectancy is about 2 years, keep your equipment and accessories free of dust and dirt, and follow these rules:
- A. Please clean the equipment before use according to chapter 6.2;Remove the batteries inside the battery cassette if the equipment will not be operated for a long time;
  - B. Replace the batteries in time when the battery voltage indicate lamps were empty;
  - C. It is recommended that the equipment should be kept in a dry environment with no corrosive gases and good ventilation. The moisture and high-light environments will affect its lifetime and even might damage the equipment.
  - D. It is best to preserve the product in a place where the temperature is between -20 to 60°C and the relative humidity is less than 95%.
  - E. The packed equipment can be transported by ordinary conveyance. The equipment not be transported mixed with toxic, harmful, corrosive materials.

WARNING

- No modification of this equipment is allowed.

7.2 Disposal

Dispose of the pulse oximeter in accordance with local environment and waste disposal laws and regulations.

8 Cleaning / Disinfection

CAUTIONS

- Never immerse or soak the oximeter.
- We recommend that the oximeter be cleaned and disinfected after use every time or determined by your hospital's policy, to avoid long term damage to the oximeter.
- Never use cleaning agents/disinfectants other than the recommended.
- The sensor component is not cleaned and disinfected during testing.

8.1 Cleaning

- The recommended cleaning agents include: water
- A. Shut down the pulse oximeter and remove the battery.
  - B. Clean the oximeter with cotton or soft cloth moistened with water.
  - C. After cleaning, wipe off the water with a soft cloth.
  - D. Allow the oximeter to air dry.

8.2 Disinfection

- The recommended disinfectants include: ethanol 70%, isopropanol 70%.
- A. Shut down the pulse oximeter and remove the battery.
  - B. Clean the oximeter as instructed above.
  - C. Disinfect the oximeter with cotton or soft cloth moistened with one of the recommended disinfectants.
  - D. After disinfection, be sure to wipe off the disinfectant left on the oximeter with a soft cloth moistened with water.
  - E. Allow the oximeter to air dry.

9 Accessories

- One lanyard.
- Two AAA batteries.
- One user manual
- One certificate card

Note:

- For particular configuration of accessories please refer to the product package list.

10 Troubleshooting

10.1 Troubleshooting

WARNINGS


- Necessary maintenance must be performed by qualified service personal ONLY.
- Users are NOT permitted to maintain the equipment by themselves.
- There are NO replaceable components in the equipment.

Trouble	Possible Reason	Solution
The oximeter can't turn to measure mode	The batteries are completely exhausted	Please replace the batteries
	An incorrect battery installation.	Verify and correct the batteries installation
	The oximeter breakdown	Please contact local service
The display is off suddenly.	The device will turn into sleep mode automatically if there is no signal in 8 seconds	Press the button again to reactivate the oximeter
	The batteries are completely exhausted	Please replace batteries
The SpO2 and Pulse Rate display unstable	The luminescent or photoelectric window is sheltered by some object.	Check the luminescent and photoelectric window.
	Excessive movement	Stop moving finger, hand and body
	The finger is not placed inside deep enough.	Place the finger properly and try again.
	Finger size is not within the recommended range	Change another finger
	Excessive ambient light	Avoid the excessive light
	Pulse rate value of the cyclical fluctuations.	The measurement is normal, and the patient is arrhythmia.
The SpO2 and PR are not displayed normally	The finger is not properly positioned	Place the finger properly and try again.
	The patient's SpO2 is too low to be detected	Try again, Go to a hospital for a diagnosis if you are sure the device works all right

Appendix A

The equipment complies with the requirement of standard EN 60601-1-2 "Electromagnetic Compatibility - Medical Electrical Equipment"

Guidance and manufacturer's declaration – electromagnetic immunity			
The MD1913 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the MD1913 Pulse Oximeter should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment –guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±8kV, ±15kV air	±8kV contact ±8kV, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTEUT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The MD1913 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the MD1913 Pulse Oximeter should ensure that it is used in such an environment.			
Immunitytest	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the MD1913 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10V/m	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80MHz to 800 MHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800MHz to 2.7 GHz  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnet 
NOTE 1 At 80 MHz and 800 MHz, the higher frequencyrange applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration – electromagnetic emissions		
The MD1913 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the MD1913 Pulse Oximeter should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The MD1913 Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MD1913 Pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations flicker emissions IEC 61000-3-3	N/A	

Recommended separation distances between portable and mobile RF communications equipment and the Medical MD1913 PULSE OXIMETER			
The MD1913 Pulse Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MD1913 Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MD1913 Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0,01	0,12	0,04	0,07
0,1	0,37	0,12	0,23
1	1,17	0,35	0,7
10	3,7	1,11	2,22
100	11,7	3,5	7,0
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.  NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			