



# V1 Compressor Nebuliser

Effective treatment for respiratory conditions



## FACT BOX

**Product Code:**  
MD2006

**Product Category:**  
Diagnostics.

**Product Approvals:**  
CE marked.

**Accessories Included:**  
x1 Atomiser  
x1 Mouthpiece  
x1 Connecting Tube (1.5m)  
x1 Adult Mask  
x1 Child Mask  
x1 Nosepiece  
x10 Filters



The Medicare V1 Compressor Nebuliser is a robust nebuliser ideal for home based, frequent use by adults or children. This highly effective nebuliser combines a short inhalation time with intelligent medication delivery thanks to its powerful compressor. The compressor generates a steam of compressed air which breaks up liquid medication and converts it into a mist so it can be easily and more effectively inhaled by the patient through the mask, nosepiece or mouthpiece.

## FEATURES

- Easy to use.
- Short inhalation time.
- Ideal for home based, frequent use by adults and children.

## DIRECTIONS:

1. Hold the medicine cup upright and turn on the nebuliser.
2. Inhale the medication as instructed by your physician.
3. Turn off the nebuliser once you are finished with the medication.

## RECOMMENDED USE:

Suitable for people with respiratory illnesses like Asthma, COPD & Cystic Fibrosis.