



# Medicare Quick-Release Tourniquet

**Product Code:**  
MD1840

## DESCRIPTION:

- Quick Release Elasticated Tourniquet for the control of bleeding.
- Controls bleeding by stopping the flow of blood to the wound.
- One hand operation allows no pinching, jamming or catching.
- Quick release buckle provides instant removal.

## CLEANING:

- To clean, wipe with a damp cloth or wash in soap and water.
- Leave to air-dry.

## STERILISATION:

- A soft hand brush can be used to remove encrusted material.
- The cuff may also be wiped with isopropyl alcohol. All tubing should be cleaned, rinsed, and dried between patients and before storage, using a detergent/intermediate-level disinfectant.

