

## Ultimate Jumpers Knee Support

Improves patellar tracking and elevation by applying mild pressure on the tendon below the kneecap



### ULTIMATE



The tubular stay allows for uniform pressure to help ease pain and discomfort.



Adjustable tension strap provides controlled compression.



Lined with breathable fabric which moves moisture away from your body keeping you cool and comfortable.



Product Code:

**PX304**

Composition:

65% AIRFLEX (Neoprene), 15% Polyester, 20% Nylon.

Recognised by



Official Partner to

