



SCHILLER
The Art of Diagnostics

SCHILLER BR-102 PLUS AMBULATORY BP RECORDER

The Schiller BR-102 Plus Ambulatory BP Recorder is a discreet, compact, long term blood pressure monitoring unit which provides the highest level of accuracy yet lets the patient go about their daily lives undisturbed. Using both auscultatory and oscillometric measurement, the motion-tolerant technology prevents re-inflation and failed readings.

FEATURES:

- Recording over 24 or 48 hours.
- Voice recording of patient data.
- Examination of adults and children thanks to a special algorithm.
- Oscillometric & auscultatory measurement.
- Cuff only inflates as much as needed, the unit can barely be seen, felt or heard during use.
- Free programming and storage up to 4 individual measurement programmes.

SOFTWARE FEATURES:

- Patient diary with automatically entered patient data.
- Integration with the SCHILLER Data Management System SEMA-200.
- E-mail transmission.
- EMR system interface (GDT).
- PDF generation.
- Serial comparison of different recordings.



Fleming Medical can offer Capital Leasing, Servicing, Training and additional accessories with this product, please contact us to request a quote.

PRODUCT DETAILS

CODE	DESCRIPTION
SHL430	SCHILLER BR-102 PLUS AMBULATORY BP RECORDER