



Blood Pressure Monitors

LIFESense®

A1 Blood Pressure Monitor	A2 Plus Blood Pressure Monitor	A5 Blood Pressure Monitor	A7 Plus Blood Pressure Monitor
MD1814	MD1803	MD1806	MD1807
Non Bluetooth Device - manually enter results on Medicare LifeSense® App.	Bluetooth Enabled - sync results to Medicare LifeSense® App on your smart device.	Non Bluetooth Device - manually enter results on Medicare LifeSense® App.	Bluetooth Enabled - sync results to Medicare LifeSense® App on your smart device.
One touch measurement.	One touch measurement.	One touch measurement.	One touch measurement.
Irregular Pulse Rate Detection.	Irregular Pulse Rate Detection.	Irregular Pulse Rate Detection.	Irregular Pulse Rate Detection.
WHO Risk Category Indicator.	WHO Risk Category Indicator.	WHO Risk Category Indicator.	WHO Risk Category Indicator.
Users			
2 x users, 60 memories per user.	2 x users, 99 memories per user.	2 x users, 60 memories per user.	2 x users, 250 memory groups per user and guest option.
Special Features			
Digital LCD display - (60mm x 41mm).	Large easy-to-read display - (60mm x 77mm).	Extra Large Full Colour Screen - (120mm x 80mm).	Extended LCD Display - (124mm x 76mm).
Large soft touch buttons and text for high visibility.	Soft touch buttons and Large text for high visibility.	Large buttons and text for high visibility.	Automatic Cuff Wrap Detection.
		Automatically averages last 3 measurements.	Blood Pressure Trend Graph - Daily/Weekly BP Results.
			On Screen Pulse Rate Display.
			Audible results option.
Cuff Circumference			
Cuff: 22cm - 32cm.	Cuff: 22cm - 32cm.	Large cuff size: 22-42cm.	Large cuff size: 22-45cm.
Power and Adaptor			
4 x AAA batteries (incl.)	4 x AA batteries (incl.)	4 x AA batteries (incl.)	Charging Port (Type-C).
AC adaptor and USB cable (not included).	USB 5V cable (optional).	AC adaptor (not included).	



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid

Ref Code: MDPOS304