



Maternity Pads

Absorbent pads with hydrophobic backing



FACT BOX
Product Code:
MD4861

Product Category:
Home Healthcare

Contents:

Maternity Pad is composed of 4 layers:

- 100% pure cellulose fluff and absorbent layer (bleached without the use of chlorine)
- Layer of cellulose paper 19gm/m²
- 2 white sheets of non-woven 18gm/m² paper
- non-woven hydrophobic paper 28gm/m² sealed with hotmelt

Medicare Non-Sterile Maternity Pads are specially designed towels for new mothers. These soft and absorbent thick maternity pads cushion and protect, offering maximum comfort while featuring a unique slim-line design.

FEATURES:

- 10 pads per pack.
- Ultra absorbent.
- Slim design.
- No lateral leaks.

DIRECTIONS:

Ideal for post birth use.

Perfect for night time use and to help keep your skin feeling dry and healthy.

CONTRAINDICATIONS AND PRECAUTIONS:

For hygiene reasons and to promote healing, change pads frequently.