



AED Wall Cabinet with Alarm and Strobe Light



FACT BOX

Product Code:
WAP-812-A01

Product Category:
Defibrillation

Material:
Aluminium

Product Approvals:
CE Approved

Our AED Wall Cabinet comes equipped with both an audible and visual Alarm System.

Designed to help prevent AED theft and alert people nearby to an emergency situation.

When the cabinet door is opened, a loud audible alarm and flashing light are activated to draw immediate attention.

FEATURES:

- Material: Aluminium.
- Size: H:38cm x W:35cm x D:18.5cm.
- Features Alarm and Strobe Light.
- Waterproof level: IP 55.
- 2 years warranty.
- Weight: 3.4kgs.
- Ventilation and Heat Dissipation System: Built-in ventilation design helps regulate internal temperature.
- IP55 Waterproof & Dustproof Protection: Rated IP55, providing reliable protection for your AED.
- Rust and UV-Resistant Aluminum Body.
- Eco-Friendly Coating & Printing.

INDICATIONS FOR USE:

Lightweight and Easy to Install: Aluminum structure ensures easy installation and maintenance for storage and protection of suitable AED Defibrillators.



WWW.FLEMINGMEDICAL.IE

Freephone: 1800 307777 Email: healthcare@flemingmedical.ie