



Outdoor AED Wall Cabinet with Heater and Keypad Lock



FACT BOX

Product Code:
WAP-812-M3KH

Product Category:
Defibrillation

Material:
Metal 2.0mm+PC glass

Product Approvals:
CE Approved

Our Outdoor AED Wall Cabinet comes equipped with a mechanical combination lock.

Provides secure access while remaining simple and reliable during emergencies.

FEATURES:

- Material: Metal 2.0mm+PC glass.
- Size: 48cm x 38cm x 24cm.
- Waterproof level: IP 56.
- Features Alarm, Heater and Keypad Lock.
- 2 years warranty.
- Weight: 14.2kgs.
- Built-in Heater with Temperature Control – Maintains optimal internal temperature, ensuring reliable AED performance even in cold conditions.
- IP56 Waterproof & Dustproof Protection: Weather-resistant construction ensuring reliable performance in harsh outdoor environments.
- Mechanical Combination Lock – Provides secure access while remaining simple and reliable during emergencies.
- Eco-Friendly Coating & Printing.

INDICATIONS FOR USE:

Lightweight and Easy to Install: Metal structure ensures easy installation and maintenance for storage and protection of suitable AED Defibrillators.



WWW.FLEMINGMEDICAL.IE

Freephone: 1800 307777 Email: healthcare@flemingmedical.ie