



Peak Flow Mouthpieces

Essential hygiene support for peak flow monitoring.



FACT BOX

Product Code:
MD1012

Product Category:
Diagnostics - Respiratory Health

Contents:
50 x Disposable Mouthpieces

Caution:
Single-patient, single-use only.
Do not clean, wash, or sterilise.
Inspect before use, discard if damaged or deformed.
Store in a clean, dry environment.

LifeSense® Breathe Peak Flow Mouthpieces are designed to support hygienic, accurate respiratory assessment in pharmacy consultation rooms and clinical environments.

Ideal for pharmacies offering asthma reviews, COPD monitoring, or routine respiratory care, this pack of 50 provides a convenient, cost-effective consumable solution for professional use.

FEATURES:

- Single-use design for optimal hygiene in consultation rooms
- Compatible with the LifeSense® Breathe Peak Flow Meter
- Made from PP5 (polypropylene) for safe, durable, hygienic use
- Reduces cross-contamination risk during respiratory assessment

DIRECTIONS:

1. Remove one disposable mouthpiece from the pack.
2. Attach the mouthpiece securely to the LifeSense® Breathe Peak Flow Meter.
3. Perform peak flow measurement as directed in the device instructions.
4. Remove and dispose of the mouthpiece after a single test or single patient interaction.

CONTRAINDICATIONS AND PRECAUTIONS:

These mouthpieces are single-use only and intended to maintain hygienic conditions during peak expiratory flow testing. Do not reuse or share mouthpieces between patients. Dispose of immediately after use in accordance with local clinical waste guidelines. Ensure the mouthpiece is securely fitted to the compatible peak flow meter before testing.