



# ADVANCING HEALTHCARE SOLUTIONS



## HEALTHCARE PRODUCT CATALOGUE

[FLEMINGMEDICAL.IE](http://FLEMINGMEDICAL.IE)



# Welcome to Fleming Medical

## Your Trusted Partner

Fleming Medical has been a trusted partner to Irish healthcare professionals for over 35 years providing the highest quality medical consumables, specialist devices and lifesaving equipment.

We believe that it is no longer sufficient to simply sell healthcare products. Understanding the needs of the healthcare professional and patient, providing a hands-on after-sales service, and finding the most effective solution is essential to our business.


## Our Customer Commitment



Customer Service  
from Product Experts



Equipment Servicing  
Team & Calibration



24 hour delivery  
via DPD



Certified to ISO 13485:2016  
by TUV Rheinland



In-House Quality  
& Regulatory Team



CPD Training & Education  
On-Site & Online

## World-Leading Brands We Work With





# FLEMING MEDICAL ACADEMY

## Professional Skills Training Platform

### Online Training Courses

Our engaging Continuing Professional Development (CPD) virtual training modules include First Aid and an interactive Basic Life Support (BLS) simulation to give a “hands-on” feeling to the online course.

#### Course Features:



Interactive  
& Self-Led



Ideal for professionals  
with busy schedules



12 Month Access to  
Course materials



Certificate of Completion  
provided when completed  
successfully



### Healthcare Product Training




#### Powerheart G5 Defibrillator

- Powerheart G5 Key Features
- How to safely operate the Powerheart G5
- Correct application of the electrodes
- Recommended CPR Technique



#### BR102 24hr Blood Pressure Monitor

- BR102 Key Features
- Correct application of the cuff
- Operation of the device
- How to use the software



#### FRED PA-1 Defibrillator

- FRED PA-1 Key Features
- How to safely operate the FRED PA-1
- Correct application of the electrodes
- Recommended CPR Technique



Scan the QR Code here to link to  
Training Platform at [www.flemingmedical.ie/training](http://www.flemingmedical.ie/training)



# CONTENT

## DEFIBRILLATION

Schiller FRED Easyport Defibrillator	9	Schiller FRED Easyport Accessories	15
Schiller FRED PA-1 AED and Trainer	10 - 11	Schiller FRED PA-1 AED Accessories	16
Schiller Defigard Touch 7	12	Cardiac Science G5 AED Accessories	17
Cardiac Science G5 AED and Trainer	13 - 14	Wall Cabinets and Bracket	18

## DIAGNOSTICS - SPECIALIST DEVICES

Schiller BR-102 Plus BPM	20	SECA MBCA 525 Scales	30
Schiller Medilog AR Holter Recorder	21	Patient Measurements	31-32
Schiller Cardiovit AT-102	22	SECA Baby & Toddler Scales	33
Schiller Cardiovit FT-1	23	SECA Rucksack	33
Schiller Cardiovit CS-104	24	SECA Baby Measuring Mat	34
Schiller Cardiovit CS-200 Excellence	25	SECA Roll Up Measuring Tape	34
Schiller Cardiovit AT-180	26	Lasso Tape	34
Aquarius Pro Patient Monitor	27	Schiller BR-102 Plus BPM Accessories	35-36
Aquarius Advanced Patient Monitor	27	Aquarius Patient Monitor Accessories	37
EDAN SE-301 ECG	28	EDAN SE-301 ECG Accessories	38
Schiller Ganshorn Spiroscout	29		

## GENERAL DIAGNOSTICS

TIMESCO Sphygmomanometer	40	Medicare LifeSense® Universal BPM Cuff	45
ERKA Sphygmomanometer	40	Medicare Compressor Nebuliser	46
TIMESCO Otoscope	41	Medicare Infrared Ear Thermometer	46
TIMESCO Ophthalmoscope	41	Medicare Prove Covers	46
TIMESCO Orion Standard Deluxe Set	41	Medicare Non-Contact Thermometer	47
TIMESCO Stethoscopes	42	Medicare LifeSense® T2 Thermometer	47
SECA Floor & Digital Scales	43	EDAN SD1 Foetal Pocket Doppler	48
Medicare LifeSense® A5 BPM	44	EDAN Pulse Oximeter	48
Medicare LifeSense® A4 BPM	45	Medicare LifeSense® Pulse Oximeter	48

## FURNITURE

Sunflower Gas Lift Stool	50	Sunflower Records Trolley	52
Sunflower Single Tier Couch Step	50	Sunflower Surgical Trolleys	52
Colenso Couch Upholstery	50	Sunflower Four Panel Screen	52
Sunflower Vista Trolleys	51	Telescopic IV Pole	52

## SIMULATION

Life/form® Adult CRiSis™ Manikin	54	Simulaids® Brad Adult™ CPR Manikin	56
Full-Body Trainer	54	Simulaids® Brad™ VTA	56
Adult Choking Manikin	54	Life/form® GERi™ Nursing Skills Manikin	57
Instructor's Economy Starter Package	55	Life/form® KERi™ Nursing Skills Manikin	57
EMT Casualty Simulation Kit	55	Life/form® TERi™ Geriatric Trainers	58
Tension Pneumothorax Simulator	55	Obstetrical Manikin with Overlays	59

# CONTENT

## SIMULATION

Rescue Randy Manikin	59	Demo Dose® Clear Ampules	64
Life/form® Basic Child CRiSis™ Manikin	60	Demo Dose® Simualted Red Blood Cells	64
Life/form® Deluxe Child CRiSis™ Manikin	60	Demo Dose® Mini Vial	65
STAT Baby Advanced	61	Demo Dose® Nitr Spray	65
Life/form® C.H.A.R.L.I.E Resus Simulator	61	Demo Dose® Code Drug II Bundle	65
Life/form® Venipuncture & Injection Arm	62	Demo Dose® Powder Vial	66
Life/form® Injection Simulator	63	Demo Dose® Inject-Ed™	66
Peritoneal Dialysis Simulator	63		

## FIRST AID

First Aid Kits	68 - 70	Sterile Saline Solution	73
CPR Face Shield	70	Emergency Eye Wash	73
Burn Gel & Dressings	70 - 61	First Aid Supplies	74
Burn Plasters	71	Skin Sanitiser Spray	74
Spray Plaster	72	First Aid Scissors	75
Alcohol Prep Pads	72	Plasters	76
Reusable Hot/Cold Pack	72		

## BANDAGES

Conforming Bandages	78	Tube Bandages	82
Elastic Adhesive Bandages	78	Stockinettes	82 - 83
Elastic Crepe Bandages	79	Net Bandages	83
Elastic Support Bandages	79	Cotton Bandages	83
Clip Fasteners	80	Tubular Bandages	84
Cohesive Bandages	80	Tubular Bandage Applicators	85
Compression Bandages	81	Foam Collar & Cuff	85
Gauze Bandages	81	Gauze Rolls	85
Muslin Bandages	81	Orthopaedic Padding	86

## ADVANCED WOUNDCARE

Adhesive & Surgical Dressings	90	Ribbon Gauze	95
Foam Dressings	90 - 91	Absorbent Pads	95 - 96
Gauze Swabs	92	Non-Adherent Dressings	96
Non-Woven Swabs	92	Sterile Eye Pads	96
Wound Pads	92	Silk Tape	97
Transparent Dressings	93	Paper Tape	97
Closure Strips	93	Transparent Tape	97
Paraffin Gauze Dressings	94	Zinc Oxide Stretch Tape	98
Absorbent Dressings	94	Zinc Oxide Tape	98
Dressing Tapes	94		

# CONTENT

## PROCEDURE PACKS

Dressing Sets	100	Suture Pack	102
Catheterisation Sets	101	Emergency Delivery Set	102
Mouthcare Set	101		

## SURGICAL

Theatre Gowns	104	X-Ray Gauze Swabs	106
Scrub Suits	104	Tracheotomy Sponges	107
Examination & Visitor Gowns	105	X-Ray Laparotomy Sponges	107
Surgical Drapes	105	Kidney Trays	107
Disposable PE Shoes	105	X-Ray Detectable Gauze Balls	108
X-Ray Abdominal Swabs	106	Sterile & Non-Sterile Gauze Balls	108
X-Ray Detectable Sterile Swabs	106		

## INFUSION

Hypodermic Needles	110	Syringe Luer Slip without Needle	112
Hypodermic Syringe	110	Catheter Syringe without Needle	112
Solution Administration Set	110	Two-Way IV Cannula	112
Syringe Luer Lock without Needle	111		

## OXYGEN MASKS & AIRWAYS

Oxygen Mask with Tubing	114	Suction Tubing	115
Disposable Nebuliser Kit	114	Nasal Oxygen Cannula	116
Non-Breather Mask	115	Guedel Airways	116
Yankauer Handle	115		

## PATIENT CARE

Cotton Wool Roll	118	Disposable Slippers	119
Cotton Balls	118	Sick Bags	120
Non-Woven Dry Wipes	118	Sick Bag Dispenser	120
Disposable & Impermeable Bibs	119	Hands & Face Personal Washcloths	120

## MATERNAL CARE

Maternity Pads	122	Sterile & Non-Sterile Baby Blankets	122
Paediatric Urine Collector	122		

## EMERGENCY

Evacuation Sheet	124	Suction Liners	126
Spider Straps	124	BVM Resuscitator	126
Head Immobiliser Block	124	Community Care Bag	126
Suction Units	125 - 126		

# CONTENT

## INFECTION CONTROL

Latex Powdered Gloves	128	Disposable PE Overshoes	130
Latex Powder-Free Gloves	128	Polyethene Aprons	130
Nitrile Powder-Free Gloves	128	Protect Barrier Gown	131
Glove Dispenser Rack	128	Cleaning Cloths	131
Hand Sanitiser	129	Face Shield	131
Detergent & Disinfectant Wipes	129	Safety Glasses	132
Disposable Face Masks	130	HSE Metal Bins	132

## MEDICINE STORAGE

Fridges	134	Zip Lock Bags	135
Hygrometer	134	Pill Boxes	136 - 138
Sample Containers	135	Pill Splitters/Crushers	138

## CONTINENCE CARE

Attends® Pull-Ons	140	MoliCare® Premium Pull-Up Pants	142
Attends® Cover-Dri Underpads	141	MoliCare® Premium Bed Mats	142

## MEDICAL CONSUMABLES

Disposable Pen Torches	144	Drawsheets	145
Gallipots	144	Examination Couch Covers	146
Tongue Depressors	144	Examination Table Covers	146
Foam Cleaning Sticks	145	Mattress Covers	146
Crepe Paper	145	Heavy-Weight Body Bag	146



DEFIBRILLATION

## BE RESCUE READY

We work with world-leading automated external defibrillator (AEDs) brands, Schiller, Zoll and Cardiac Science, who are trusted by first responders, healthcare sites, businesses and communities to treat sudden cardiac arrest (SCA). Our AEDs are designed with everyday rescuers in mind.

World-Leading AED Brands We Work With



## After-Sales AED Programme

As suppliers of vital and sometimes life-saving medical equipment, we know that excellent product training, servicing, and calibration are essential to ensure that these devices work correctly when they are needed and are operated properly. That's why we don't just supply medical equipment and devices, we also provide a robust after-sales programme to ensure the operational precision of all of our AEDs.

Our Medical Equipment Service Team calibrates, services and maintains life-saving and medical diagnostic equipment for clinical, private and public environments.



## Benefits

- ✓ Free Product Demos.
- ✓ Full Range of Accessories & Spare Parts.
- ✓ Online & In-Person AED Training.
- ✓ Manufacturer Trained Repairs.
- ✓ Maintenance Contracts Available.

## In-House Service Department

Direct contact with our Equipment Service Team is available to our customers.







## Schiller Fred Easypart Plus Automatic Portable Defibrillator

The Schiller Fred Easypart Plus Automatic is a small, lightweight, portable defibrillator device. At-risk patients can carry their own rescue device after they and their families have been trained by their doctor. This dramatically reduces the response time for treating ventricular fibrillation and tachycardias, offering victims a much better chance of survival.

- Biphasic pulsed defibrillation impulse (Multipulse Biowave®) - pulsed energy which is more gentle on heart tissue.
- High resolution LCD screen.
- Battery will last up to 3 years in standby mode\*.
- Charging time from shock recommended to shock standby less than 10 seconds.
- Extremely short shock recovery.
- Weighs 490g including batteries.
- Dimensions: 133 x 126 x 35mm.

### CODE

SHL800

NOTE: 0% VAT on all AEDs



Lightweight  
& compact



Training materials  
available



Long  
battery life

\*We recommend to change it when it reaches 3 years.





### Schiller Fred PA-1 AED

The Schiller FRED PA-1 Defibrillator starts up immediately and guides the rescuer step by step during the entire resuscitation process. It is embedded with the FreeCPR feature, the only system that helps the rescuer to reach the proper compression frequency without any additional hardware.

- Analyses the patient's heart rhythm and determines whether a shock is necessary.
- Prompts the rescuer to deliver the shock by pressing the shock button.
- FRED PA-1 immediately advises to call the emergency number directly listed on the device.
- Chest compression assistant:
  - Metronome.
  - FreeCPR: Feedback on the frequency of the chest compression.
- Two types of electrodes, for adults and children.
- 10 Year Warranty.

#### CODE

SHL776

NOTE: 0% VAT on all AEDs



#### Fast Set-up

Pre-connected electrodes allow faster application



#### Semi-Automatic

The FRED PA-1 prompts the rescuer to deliver the shock by pressing the shock button



#### Battery Shelf Life

Medical grade 6 year battery



**SCHILLER**  
The Art of Diagnostics



FREE  
Demonstration  
Available

## Schiller FRED PA-1 AED Trainer

A training device that can realistically train the user of an AED. Several scenarios make it possible to simulate different resuscitation situations. Associated to a manikin, the FRED PA-1 Trainer provides training for the whole resuscitation process.

- Repositionable training electrodes to train the positioning of the electrodes on the torso.
- The FRED PA-1 Trainer delivers detailed step-by-step instructions to the user.
- During the chest compression phase performed on the manikin, the FRED PA-1 Trainer indicates the rhythm to follow using the metronome.
- Several languages available.
- The Schiller training app is available in the Apple App Store and the Google Play Store.

### CODE

SHL783



**SCHILLER**  
The Art of Diagnostics

NOTE: 0% VAT on all AEDs





### Schiller Defigard Touch 7

The Schiller Defigard Touch 7 is an indispensable tool for rescuers. The device is extremely compact and offers the latest defibrillation technology in combination with comprehensive monitoring functions. It is the first emergency monitor/defibrillator equipped with a touch screen, making it the most intuitive device on the market, as well as with the latest data transmission technology.

- 7" high definition touch screen device.
- Touch'n'Save user interface.
- ECG with 6 to 12 leads.
- SpO<sub>2</sub>, NIBP, IBP, EtCO<sub>2</sub>.
- AED, Manual Defibrillator and temperature.
- Connectivity via cellular network and WiFi.
- Working on a built in Li-ion battery.

#### CODE

SHL753

NOTE: 0% VAT on all AEDs



High definition  
touch screen



Ambulance  
compatibility



Can print wireless on  
a thermal printer





## Cardiac Science Powerheart G5 AED

Whether you are an experienced rescuer or a first-time responder, the Powerheart G5 automated external defibrillator (AED) provides a powerful combination of features that help rescuers provide sudden cardiac arrest (SCA) victims with swift, effective, and life-saving therapy.

- First AED to combine variable escalating energy for customised patient therapy.
- Faster shock times and real time CPR feedback.
- Durable design tested to military standards with an IP55 rating protection from dust and water.
- Highly visible "Rescue Ready" indicator confirms rigorous daily, weekly and monthly self-tests.
- Provides customised shocks tailored to each patients needs.
- 8 year warranty.

CODE

CSG502

NOTE: 0% VAT on all AEDs



User-paced instruction  
from RescueCoach™



Indicator communicates  
rescue status



Intellisense™  
CPR Feedback







### Cardiac Science Powerheart G5 AED - Defibrillator Trainer Kit

The Powerheart G5 AED Trainer - the best way to prepare is realistic, hands-on defibrillator rescue experience. The Powerheart G5 AED Trainer is designed to closely mimic the Powerheart G5 AED for effective, hands-on defib training.

Each AED Trainer package includes:

- (1) AED Trainer.
- (1) Set of Adult Training Pads.
- (1) Wireless Remote.
- (1) Language CD with User Manuals.

CODE

CSG503

NOTE: 0% VAT on all AEDs







**SCHILLER**  
The Art of Diagnostics

## Schiller Fred Easyport Portable Defibrillator Accessories



**Schiller Adult Pads**  
SHL07/AP



**Schiller Paediatric Pads**  
SHL787



**Schiller Battery**  
SHL07/B

NOTE: 0% VAT on all AEDs Parts & Accessories



**SCHILLER**  
The Art of Diagnostics

## Schiller Fred PA-1 AED Accessories



Schiller Pre-Connected Pads\*  
SHL752/1



Schiller Paediatric Pads\*  
SHL787



Schiller Non-Rechargeable  
Lithium Battery SHL78



Defigard Battery SHL900

NOTE: 0% VAT on all AEDs Parts & Accessories

\*Compatible with Schiller Fred PA-1 and Schiller Defigard Touch 7.



## Cardiac Science Powerheart G5 AED Accessories



**Adult Defib Pads**  
CSG504



**Adult Defib Pads**  
CSG505



**Paediatric Defib Pads**  
CSG506



**Adult Defib Training Pads**  
CSG507



**Adult Defib Training Pads**  
CSG508



**4 Year Battery**  
CSG510



**Premium Carry Case**  
CSG511



**Ready Kit**  
CSG512



**Wall Storage Sleeve**  
FR9159

NOTE: 0% VAT on all AEDs Parts & Accessories



**Defibrillator 1000  
Mild Steel Cabinet**  
CS216

- Suitable for coastal locations, protected against corrosion with an outdoor yellow powder coat.
- Internal LED magnetic light that activates in darkness on a motion sensor and can be used as a removable torch.



**AIVIA AED External Wall  
Cabinet with Keypad**  
CS210

- The front face of the cabinet shows the survival sequence and 3 clear pictograms with instructions for use.
- Day/night lighting of the AED compartment by LED with automatic luminosity detection.



**AIVIA AED Wall Mount  
Cabinet**  
CS212

- The front face of the cabinet shows the survival sequence and 3 clear pictograms with instructions for use.
- Day/night lighting of the AED compartment by LED with automatic luminosity detection.



**Cardiac Science  
Powerheart AED Cabinet  
& Strobe Light Alarm  
Security Enable**  
FR50-00392-30

- Powerheart wall mount storage case with strobe light alarm is a highly visible, well-marked cabinet to house your facility's Powerheart AED.



**Schiller Wall Bracket  
for Fred PA-1**  
SHL34

**AIVIA**



NOTE: 0% VAT on all AEDs Parts & Accessories





## DIAGNOSTICS - SPECIALIST DEVICES

At Fleming Medical we offer our customers the latest and most innovative specialist diagnostics devices. We are proud to be a leading distributor partner with world-leading brands including Schiller, Edan and Aquarius brands.

We supply ECGs and accessories, patient monitors, and paediatric health check equipment. Our specialist devices are backed by excellent customer service and after-sales care from our service team.

## 24HR AMBULATORY BP MONITOR

The Schiller BR-102 Plus Ambulatory BP Recorder is a discreet, compact, long term blood pressure monitoring unit which provides the highest level of accuracy while allowing the patient to go about their daily lives undisturbed.

Using both auscultatory and oscillometric measurement, the motion-tolerant technology prevents re-inflation and failed readings.



Ranked on the HSE  
Dynamic Purchasing  
System for Non-Invasive  
Cardiology Equipment

### Schiller BR-102 Plus Ambulatory BP Recorder

- Recording for 24 or 48 hours.
- Voice recording of patient data.
- Examination of adults and children thanks to a special algorithm.
- Oscillometric and auscultatory measurement.
- Cuff only inflates as much as needed, the unit can barely be seen, felt or heard during use.
- Free programming and storage up to 4 individual measurement programmes.
- On current HSE Contract HSE6909.

#### SOFTWARE FEATURES:

- Patient diary with automatically entered patient data.
- Integration with the SCHILLER Data Management System SEMA-200.
- E-mail transmission.
- EMR system interface (GDT).
- PDF generation.
- Serial comparison of different recordings.



Lightweight, compact  
& discreet



Built-in pump, making  
it almost silent



Easy menu guidance with  
two control buttons

CODE

SHL430



**SCHILLER**  
The Art of Diagnostics



## HOLTER RECORDER

## SPECIALIST DEVICES



**RANKED 1ST** on the  
HSE Dynamic Purchasing  
System for Non-Invasive  
Cardiology Equipment

### Schiller Medilog®AR Holter Recorder

The Schiller Medilog®AR Holter Recorder is a robust, splash and shockproof holter that is easy to clean. Thanks to the dual-battery concept, the patients can be screened for more than 14 days without changing batteries.

- Zero-second atrial fibrillation atrial flutter detection based on true P-wave analysis.
- Analysis takes less than 90 seconds with DARWIN2 software.
- Capable of open-ended recording with dual battery function.
- 32,000Hz sampling rate on 3 channels, for true p-wave analysis, sophisticated artefact suppression and motion detection.
- The optional SpO2 sensor, connected via Bluetooth, provides additional respiratory information.
- Ranked 1st in HSE 6909 Framework.
- Microphone to record patient ID.
- Runs on internal rechargeable battery, or normal AAA batteries if desired.
- Includes robust carrying case with various carrying options.



Robust carrying case with  
various carrying options



Microphone to record  
patient ID



Runs on an internal  
rechargeable battery, or on  
normal AAA batteries if desired

#### CODE

SHL10111





Ranked on the HSE  
Dynamic Purchasing  
System for Non-Invasive  
Cardiology Equipment

## Schiller Cardiovit AT-102 G2

The new generation 12 Lead ECG for hospitals and clinics provides superior workflow efficiency with Wi-Fi and advanced interpretation as standard.

- 12 Lead ECG.
- High resolution 8" colour display.
- Wi-Fi and LAN connectivity as standard.
- Smart battery technology guarantees over 8 hours of ECG recording.
- Unlimited product training.
- Trolley available.
- Continuous remote and on-site tech and maintenance support.
- Stand alone and CIS compatible ECG machine.

### CODE

SHL145/1



Internal high-resolution  
A4 thermal printer



Touch function keys



Sealed alphanumeric  
keyboard



**SCHILLER**  
The Art of Diagnostics



## Schiller Cardiovit FT-1

Schiller's latest innovation is ultra portable, has an intuitive 8-inch high-resolution multi-touch screen, an integrated thermal printer and is FDA approved.

Designed for users who value state-of-the-art technology, the Cardiovit FT-1 is a networked device providing connectivity and a paperless workflow.

- Power and flexibility of a PC in a portable ECG device.
- Fast and secure bidirectional Wi-Fi communication.
- Various optional diagnostic analysis programs: automatic interpretation ETM, ETM Sport based on Seattle Criteria, CCAA for early STEMI detection.
- Resting rhythm for up to 4 minutes.
- Hook-up adviser with colourcoded waveforms and anatomical model.
- Lead reversal detection.
- Outstanding signal quality for adult and paediatric ECGs.
- Worklist (optional).

CODE

SHL352



Intuitive multi-touch screen



Connectivity and a paperless workflow



Anatomic hook-up adviser



Lead reversal detection





### Schiller Cardiovit CS-104

Cardiovit CS-104 combines all essential ECG features and desired possibilities in an integrated system. Cardiovit CS-104 is also available as resting ECG PC software and spirometry PC software.

With the Cardiovit CS-104, Schiller offers a complete ECG system for practices, clinics and hospitals. The unit includes a PC, a monitor, a trolley and a wireless (MS-12 blue) or wired (MS-12 USB) acquisition module. Optionally, Cardiovit FT-1 can be used for ECG acquisition. Everything is pre-installed and tested by Schiller, ready to be used in the clinical practice.

- Resting ECG with the Schiller ECG analysis program ETM as well as ETM Sport, the automatic interpretation for athletes.
- Resting rhythm recording for up to 20 minutes, with optional arrhythmia detection.
- ECG Framer: 10-second resting ECG from resting rhythm.
- Exercise ECG with optional arrhythmia detection.
- Vector ECG, optional.
- Spirometry, based on the ultrasound technology of the SpiroScout SP plus, optional.
- Integrated workstation for data review from all Schiller devices connected to the Schiller Server.
- Bi-directional connectivity (including worklist and PDQ) to EMR, HIS and PACS systems.

#### CODE

SHL562



Exercise ECG



User friendly



Hospital



User interface





### Schiller Cardiovit CS-200 Excellence

The CS-200 Excellence cardiac assessment system for exercise testing from Schiller puts great emphasis on performance.

With industry-leading algorithms, the CS-200 Excellence enables you to quickly and easily optimise your digital workflow with seamless connectivity to EMR, SEMA cardiology information system and PACS.

- Outstanding data quality with algorithms that provide the clearest, most accurate test results.
- ECG measurement software with accurate measurements of recorded cardiac signals.
- Resting rhythm.
- Hook-up adviser with signal quality check, helping to improve poor waveform quality.
- Pharmacological exercise testing.
- 16-lead acquisition and analysis for resting and exercise ECG: based on real measurement data, not just calculations.
- Various diagnostic analysis programs: ETM, ETM Sport, HyperQ (Resting and Stress), SAECEG.

CODE

SHL320



Seamless connectivity to  
EMR, SEMA & PACS



Various diagnostic  
analysis programs



16-lead ECG acquisition  
and analysis



Pharmacological  
exercise testing



**SCHILLER**  
The Art of Diagnostics





### Schiller Cardiovit AT-180

This multi-use workstation combines the largest touch screen with advanced clinical tools.

Schiller's high-end electrocardiograph AT-180 is based on the new Schiller software platform.

It offers robust design for the busiest hospitals, as well as a combination of applications that makes it the ideal multi-use workstation for private practices.

- 18.5" high-resolution multi-touch screen for easy ECG review.
- Full-size keyboard with a durable, exchangeable cover keeping dust, dirt and liquids out, making it highly hygienic.
- Resting rhythm for up to 20 minutes.
- ECG Framer: Create a 10-seconds resting ECG from a resting rhythm recording.
- Hook-up adviser with colour-coded waveforms and anatomical model.
- Lead reversal detection.
- Fast and secure bidirectional Wi-Fi communication, PDQ
- Worklist.
- Exercise test (optional).
- Various diagnostic analysis programs: ETM, ETM Sport, SAEKG, Vector ECG, arrhythmia detection, ECG re-measurement (optional).
- Resting 16-lead acquisition and analysis: based on real measurement data, not just calculations (optional).

CODE

SHL561



High-resolution  
multi-touch screen



Fast and secure bidirectional  
Wi-Fi (communication)



Exercise ECG



Hospital



**SCHILLER**  
The Art of Diagnostics





### Aquarius Vital Signs Monitor ECG / SPO2 / NIBP / Temp

The Aquarius Patient Monitor ECG / SpO2 / NIBP / Temp is one of the most popular patient monitors due to its wireless capabilities and convenient features. Portable, lightweight and of sturdy design, it is suitable for adult, paediatric and neonatal patients. Its rechargeable Li-ion battery provides up to 4 hours of uninterrupted work.

- 5 Lead ECG with Respiratory Rate\*
- 8" Colour TFT LCD Touch Screen.
- Standard Configuration: NIBP, SpO2, ECG & Temperature.
- Flexible Parameters: Configuration for different clinical environments.
- Continuous Monitoring Mode with Alarm Function.
- Interchangeable Patient Trend Views: Graphical and Tabular.

#### OPTIONAL EXTRAS:

- Roll Stand, Carry Case.
  - FREE Demo and Training.
- \*5 Lead ECG with respiratory rate, please contact your Sales Rep to request a quote.



Portable, Lightweight,  
Sturdy Design




Suitable for Adult,  
Pediatric & Neonatal  
Patients



Upgrades Available

CODE	DESCRIPTION
NOME01	Aquarius Pro Vital Signs Monitor
NOME21	Aquarius Advanced Vital Signs Monitor

# AQUARIUS



### EDAN SE-301 3-channel ECG

The EDAN SE-301 is the new generation 3-channel ECG, which incorporates the benefits of advanced technology and innovative design in an extremely compact and lightweight device. It is specifically developed for those who require ease of portability without compromising on quality.

- High resolution colour touch screen.
- Comprehensive filters and anti-noise technology.
- Enhanced powerful technical specifications provides reliable signal quality.
- 12-lead simultaneous acquisition and display.
- Support linear/2D barcode scanner.
- Support USB flash disk and micro SD card to extend memory.
- Built-in rechargeable Li-ion battery.
- Bi-directional communication with SE-1515 Data Management System.
- Support external printer via USB port.

CODE

EDA012



Lightweight,  
Less than 1kg



Accurate &  
Reliable Results



Results are Easy to  
Print, Share & Save





### Schiller Ganshorn Spiroscout

The SpiroScout is a complete lung function laboratory based on the unique ultrasonic measurement principle of Ganshorn. It enables the exact determination of all spirometric parameters. As the first of its kind the SpiroScout simultaneously measures flow and gas density providing all the necessary information about the volume and gas exchange from one single measurement.

- High precision spirometry due to direct flow measurement.
- Realtime BTPS correction, drift free, high accuracy.
- Different options available such as rhinomanometry, oscillatory resistance, capnovolumetry, emphysema test, dead space analysis.
- No volume calibration and no warm-up time necessary.
- Perfectly hygienic concept optimising workflow and patient safety.
- No perceptible resistance and more comfort for the patients.

#### CODE

SHL356



### SECA MBCA 525 Body Composition Scales

The new compact class for mobile body composition analysis. The unique measuring mat measures all of the required parameters such as fat mass, muscle mass and body water in a fully self-sufficient way. The short cables simplify working. In addition, the mat is automatically linked to the touchscreen monitor by Wi-Fi, which can also be deactivated if required.

The large memory of the monitor has sufficient capacity for over 100,000 measurements. This enables a large number of patients and courses of therapy to be reliably documented for years. Furthermore, measurements can be reliably carried out regardless of location in hospitals or in private practices. Its lightweight, practical handle and optional mobile stand with a basket for patient files and accessories guarantees maximum mobility.

- Fast measurements executed in 17 seconds.
- Measurement of fat, water, muscle and other body parameters in lying position.
- Perfect for mobile use at only 3kg with long-lasting rechargeable battery.
- Wi-Fi connection between measuring mat and monitor guarantees reliable data transmission.
- Easy-to-understand graphic presentation of measurement results on the touchscreen monitor.

#### CODE

SEC525

**seca**<sup>®</sup>



## SECA Electronic Chair Weighing Scales

- Digital Chair Weighing Scales is well suited for hospitals and clinics, particularly geriatrics, nephrology, orthopaedics, nursing homes, retirement homes and rehabilitation centres.
- Rounded steel tube frame.
- Comfortable seat.
- Swivelling armrests - easy to move patient from bed to chair.
- Adjustable footrests.
- Locking rear wheels - safe to sit down and stand up at all times.
- Large LCD display unit behind the chair - easy to read.
- BMI facility included.
- HOLD function.
- Mains and rechargeable battery power supply.
- Weight: 24.5kg.
- Capacity 250kg/39 stone.
- Graduation: 50g <150kg > 100g.
- CE Class III Approved.
- Dimensions: 533(W) x 906(H) x 946(D) mm.

CODE

SEC955



## SECA Portable Height Measure

- At only 2.4 kg, the SECA 213 mobile stadiometer is especially suitable for use in the field but also ideal for doctors' practices and paediatric hospitals.
- The measuring rod can be dismantled into several pieces and set up easily and quickly. The spacer keeps the rod straight and stable without any fittings. As the scale is printed along the side of the measuring rod, it is easy to read off the result while measuring, thus guaranteeing precise results – up to a height of 205 cm.
- Comes with carry bag.

CODE

SEC213

**seca**<sup>®</sup>





### SECA Digital Column Scale with BMI Function

- Type: Digital.
- Capacity: 200kg.
- Graduation: 100g.
- Size: (W x H x D) 294 x 831 x 417mm.
- Weight: 6.3kg.
- Power supply: batteries or mains adapter.
- Approval: Class III.
- Functions - Tare, Hold, BMI, Weighing Range Switch Over, Auto Switch Off.

#### CODE

SEC799



### SECA Height Measure Rod for Column Scales

- Telescopic Height Measuring Rod for Column Scales.
- Measuring range: 85-200cm.
- Graduation: 1mm/ 1/8".
- Dimensions: 55 (W) x784 (H) x34(D) mm.
- Weight: 700g.

#### CODE

SEC220



### SECA Electronic Baby & Toddler Scales with Detachable Cradle

All purpose electronic baby and infant scale with detachable tray for weighing in a lying or standing position.

- Capacity 50kg.
- Graduation 20g.
- Dimensions 550 x 165 x 330mm.
- Weight 2.5kg.
- Power supply battery.

CODE
SEC385



### SECA Rucksack

- Rucksack for use with SECA Electronic Baby and Toddler Scales.

CODE
SEC427



### SECA Baby Measure Mat

- Washable and foldable.
- Fixed head piece.
- Easy to slide foot positioner.
- Measuring Range: 10-99cm/4-39".
- Graduation: 5mm/ 1/4".
- Size: 300 (W) x 140 (H) x 1340 (D) mm.

#### CODE

SEC210



### SECA Roll Up Measuring Tape with Wall Attachment

- Space-saving roll-up mechanism.
- Simple read-out of results through display window.
- Durable tape made of metal.
- The durable tape is made of metal.

#### CODE

SEC206



### SECA Baby Measuring Lasso Tape

- Hygienic disposable measuring tape of non-stretch Teflon material.
- Double function: head circumference with the front side of the tape and facial symmetry with the rear side.
- Easy-to-read, colour-coded precise results.
- Scale from 1 to 59cm.

#### CODE

SEC212

#### PACK SIZE

15



### Head Lasso-o Tape for Baby

- Head Lasso-o Tape is the recommended equipment used for measuring the circumference of a baby's head.
- Accuracy is vital and the use of the Lasso-o guarantees this.







#### CODE

HL01

#### PACK SIZE

10

### Schiller BR-102 Plus BPM Accessories

		
<p>Schiller BR-102 Plus Blood Pressure Cuffs</p>	<p>Schiller BR-102 Plus Connection Tube with Microphone SHL230/6</p>	<p>Schiller Premium Reusable Red Pouch for BR-102 Plus SHL97</p>
<p>SHL281/1 - Small - 18-26cm</p>	<p>Suited for use with Schiller BR-102 Plus small, medium and large blood pressure cuffs. Comes with microphone.</p>	
<p>SHL281 - Medium - 25-35cm</p>		
<p>SHL281/2 - Large - 35-45cm</p>		
		
<p>Schiller BR-102 Plus Charger for Ni-MH Batteries (SHL230/7) SHL230/9</p>	<p>Schiller BR-102 Plus Rechargeable Battery Ni-MH SHL230/7</p>	<p>Schiller Medilog Holter® Accessories</p> <p>Schiller Medilog AR Holter 7-Wire Cable SHL1013</p>
	<p>Charger for NiMH batteries for Schiller BR-102 Plus Ambulatory Blood Pressure Monitor.</p>	<p>Rechargeable Battery for Schiller BR-102 Plus. Ni-MH (2500mAh) battery.</p>

Schiller ECG Accessories



**Schiller Trolley X1**  
SHL96



**Schiller ECG Cable  
Support for Trolley X1**  
SHL94



**Schiller 10 Lead  
Patient Cable**  
SHL48



**Blue Sensor Adult ECG  
Electrodes** SHL132



**Bio-Adhesive Electrodes**  
SHL60



**Snap Clip Adaptor**  
SHL61



**ECG Paper** SHL107



**Battery Pack** SHL387



## Aquarius Vital Signs Monitor NOME01 and NOME21 Accessories



**Aquarius Trolley**  
NOME02



**Carry Case**  
NOME10



**Wall Mount**  
NOME01/2



**Button ECG Electrode\***  
NOME14



**Reusable Adult BPM Cuff**

Medium 27-35cm  
NOME09

Large 34-43cm  
NOME06



**Reusable Paediatric BPM Cuff**

Medium 10-19cm  
NOME08

Large 14-21.5cm  
NOME07



**Adult SpO2 Sensor**  
NOME12



**Adult Clip SpO2 Sensor**  
NOME15



**Paediatric SpO2 Sensor with 2m Long Cable**  
NOME5



**Ear Probe Covers\***  
NOME01/3



**Printer Paper\***  
NOME03



**Aquarius Rechargeable Battery** NOME17

\* FOR USE WITH THE AQUARIUS PRO VITAL SIGNS MONITOR (NOME01) ONLY.

# AQUARIUS

Edan ECG Accessories



**EDAN Rolling Stand  
(MT-201) EDA013**

Of the 5 wheels, 2 are equipped with a brake to secure the stand when it's stationary. Two baskets for the storage of medical equipment and its height can be adjusted between 58-109cm. The ECG can be oriented with a tilt of 15 degrees.



**EDAN Carry Case  
EDA014**



**EDAN SH5 Paediatric  
Sensor EDA007**

Paediatric Sensor for Pulse Oximeter. Data Management CD and RS-232 Cable.



**SpO2 Reusable Adult  
Finger Sensor (1m)  
EDA018**



**EDAN Recording Paper  
Roll EDA016**

## GENERAL DIAGNOSTICS

We partner with healthcare professionals to ensure they have access to the highest quality general diagnostic devices on the market. We provide our customers with a variety of equipment such as digital scales, pulse oximeters, stethoscopes, and blood pressure monitors with arm cuffs.

Our general diagnostic devices are backed by excellent customer service and after-sales care from our service team.

SPHYGMOMANOMETERS



TIMESCO Emerald 'Clip On' Aneroid Sphygmomanometer (2 Tubes)

- 2 tube system.
- Separate bellows and gauge tubes.
- Spring loaded clip-on rear of dial for one-handed use.
- 48mm dial for easy reading of pressure values.
- Push on tube connection (no adaptor required).
- Supplied with an adult black nylon calibrated cuff.
- Zip-up pouch for convenient storage.

CODE	SIZE
TIME200	22-32cm



ERKA Perfect Palm Sphygmomanometer

ERKA 'Perfect Palm' Sphygmomanometer is a top of the range, high quality, classic handheld blood pressure gauge.

- Ideal for use for GPs and hospitals.
- Manometer housing brass chrome-plated.
- Membrane protecting 2 tube technology.
- Display scale is 56mm diameter.
- Measuring accuracy over the total range:  $\pm 3$ mmHg.
- Screw valve for sensitive adjustment.

CODE
ERK01



## TIMESCO Orion Xenon Standard Otoscope

The Orion Xenon Otoscope has a rugged metal head with a non-reflective matt finish and provides a bright white light from a Xenon powered long-life bulb.

The otoscope head has a universal fitting for aural, nasal and also an adaptor to fit the TIMESCO Disposable Ear Speculae.

- Lifetime guarantee.

### CODE

TIME2500



## TIMESCO Orion Xenon Standard Ophthalmoscope

The Orion Xenon Ophthalmoscope features a head made from lightweight hermetically sealed durable plastics, precision optics and a latex-free rubber eyebrow rest.

- The Xenon powered bright white light bulbs, provide crystal clear illumination in ophthalmic diagnostic procedures.
- Size C 2.5v in soft pouch.
- Lifetime guarantee.

### CODE

TIME2501



## TIMESCO Orion Xenon Standard Deluxe Set

Xenon Size C 2.5v in hard case comes with:

- Interchangeable handle.
- Otoscope head.
- Ophthalmoscope head.
- Disposable speculum.
- Disposable speculum adaptor.
- Nasal speculum screw.
- Spare otoscope bulb.
- Spare ophthalmoscope bulb.
- Lifetime guarantee, including consumables.

### CODE

TIME2502





### TIMESCO Dual Head Stethoscope

- Dual head stethoscope is a solid lightweight aluminium chest piece which provides good acoustic response.
- Non-chill rim on diaphragm.
- Solid, lightweight aluminium chest piece provides good acoustic response.
- ABS ear tip for a comfortable fit.

CODE

TIME101



### TIMESCO Training Dual Head Teaching Stethoscope

- Dual head chest piece (47mm diameter).
- Designed specifically for training purposes.
- Two binaural headsets coupled with quality tubing.
- Allows a teacher and student to listen to heart and chest sounds simultaneously.
- Aluminium alloy chest piece.

CODE

TIME108



### TIMESCO Diamond Stethoscope

- Dual head chest piece.
- A high quality 48mm diameter chest piece with a floating diaphragm for excellent acoustics.
- Contains internal binaural springs.
- Supplied with a spare diaphragm and spare ear tips.
- Supplied in a foam lined box.
- Stainless steel chest piece.

CODE

TIME103



### TIMESCO Paediatric Diamond Stethoscope

- Dual head chest piece.
- Paediatric chest piece (40mm diameter).
- Infant chest piece (30mm diameter).
- Supplied with spare diaphragm and spare ear tips.
- Supplied in a foam lined box.
- Stainless steel chest piece.

CODE

TIME107



### SECA Personal Floor Scales

- Large, low level non-slip platform.
- Quick weighing.
- Reliable corrosion protection.
- Type: Mechanical.
- Weight: 4kg.
- Capacity: 150kg.
- Graduation: 1kg.
- Size: 300 (W) x 105 (H) x 460 (D) mm.
- Platform: 270 (W) x 65 (H) x 280 (D) mm.

#### CODE

SEC761



### SECA Digital Scales

- A class 3 electronic self-zeroing scales.
- Approval class: III (3).
- Capacity: 200kg.
- Graduation: 200g.
- Power supply: Battery.
- Measure: 321 (W) x 60 (H) x 356 (D) mm.
- Net weight: 4kg.
- Functions: Mobile measuring and weighing, Tap-on automatic switch-on, Automatic switch-off.

#### CODE

SEC875



### SECA Electronic Floor Scales

Scales intended for mobile use have to be lightweight. This scale is nevertheless very stable and suitable for persons weighing up to 200kg. It stands firmly on four levelling bases that guarantee stability, even when small children are weighed while held by an adult.

- Graduation Weight: 100g < 150kg > 200g.
- Weight: 4.2kg.
- Power supply: Batteries.

#### CODE

SEC877



## Medicare LifeSense® A5 Blood Pressure Monitor

Our fully automatic Medicare LifeSense® A5 Upper Arm Blood Pressure Monitor has clinically proven accuracy. The coloured screen combined with large text and buttons, makes this blood pressure monitor the ideal choice for people with poor eyesight.

- One touch measurement.
- Large buttons and text for high visibility.
- Equipped with silent pump.
- Comfortable touch interface.
- Automatically averages last 3 measurements.
- Measure during inflating technology.
- Large cuff size: 22-42cm.
- Track your blood pressure by inputting your readings manually into the Medicare LifeSense® App.

### CODE

MD1806



Irregular Heartbeat  
Detection



2 x Users  
60 Memories per user



WHO Risk  
Category Indicator  
(World Health Organisation)\*



Extra Large  
Full Colour Screen



### Medicare LifeSense® A4 Upper Arm Blood Pressure Monitor

The Medicare LifeSense® A4 Blood Pressure Monitor is designed to connect with the Medicare LifeSense® health status monitoring app so you can effortlessly monitor your blood pressure from the comfort of your own home. It uses irregular heartbeat detection and is reliable due to its clinically proven accuracy.

- Accurate measurement to  $\pm 3$ mmHg.
- Averaging of last 3 readings.
- Large cuff size: 22-42cm.
- Supersize LCD Display.
- Date and time function.
- Power: 4 x AAA batteries (included).
- USB cable (included).
- Auto shut down.



Irregular Heartbeat  
Detection



2 x Users  
60 Memories per user



WHO Risk  
Category Indicator  
(World Health Organisation)\*

CODE
MD1809



### Medicare LifeSense® Universal BPM Cuff

- Compatible with all Medicare LifeSense® Blood Pressure Monitors.

CODE	TYPE	SIZE
MD1804	Paediatric	18-26cm
MD1805	Adult	22-42cm



### Medicare V1 Compressor Nebuliser

The Medicare V1 Compressor Nebuliser is a robust nebuliser, ideal for frequent use by adults or children. This highly efficient nebuliser combines a short inhalation time with intelligent medication delivery thanks to its powerful compressor.

#### CODE

MD2006

- Suitable for people with respiratory illnesses like Asthma, COPD and Cystic Fibrosis.
- Particle size: 5µm.
- Medicine capacity: 8ml.
- Comes with full compressor nebuliser accessory kit.
- Suitable for adults and children.



### Medicare Infrared Ear Thermometer with Probe Covers

Medicare Infrared Ear Thermometer is a fast, accurate thermometer with a user-friendly ergonomic design. It uses infrared technology to measure the heat generated by the eardrum and surrounding tissue which determines the patient's temperature.

- Fast, 1 second measurement.
- Large LCD display.
- Supplied with probe covers for hygiene purposes.
- User-friendly ergonomic design.
- Fever alarm and 9 Memory spaces.
- °C/°F - Switch between Celsius and Fahrenheit.

#### CODE

MD1882



### Medicare Probe Covers for MD1882 Ear Thermometer

- 40 Medicare probe covers for MD1882 Ear Thermometer.

#### CODE

MD1883



## THERMOMETERS



**1** Fast, 1 Sec  
Reading

**3** in-1 MEASUREMENT  
Body, Surface &  
Room Temperature

### Medicare Non-Contact Infrared Thermometer

Medicare's Non-Contact Thermometer is a fast, accurate infrared thermometer with convenient features to help you monitor temperature with ease. The 3-in-1 measurement for body, surface and room temperature makes it ideal for managing temperatures around your home for babies and children.

- Fast, 1 second measurement.
- Measures without skin contact (3-5cm).
- Fever alarm.
- Large colour coded backlit display.
- 32 memory sets.
- °C/°F - Switch between Celsius and Fahrenheit.
- Silent Mode.
- Auto-off.
- Practical storage box.
- User recalibration enabled.
- Low battery indicator.

#### CODE

MD1831



### Medicare LifeSense® T2 Digital Thermometer

Medicare LifeSense® T2 Digital Thermometer for fast, accurate measuring of temperature.

- Waterproof.
- Suitable for oral, underarm or rectal use.
- 10 Second fast measurement time - 6 times faster than traditional thermometers.
- °C/°F - Celsius and Fahrenheit.
- Auto Shut-off: 10 mins after use.
- Memory recall.
- Flexible Tip.
- Fever alarm.
- 200 hr or 12 months Battery Life.

#### CODE

MD061



### Medicare LifeSense® Fingertip Pulse Oximeter

The purpose of our Medicare LifeSense® Fingertip Pulse Oximeter is to measure your SpO2 levels and pulse. This means it checks how well your heart is pumping oxygen through your body.

- Visualisation of heartbeat.
- Pulse rate sound synchronised with heartbeat.
- Alarm function (SpO2 <88%, pulse rate <40 or >120).
- One button to start & automatic power off.
- 2 AAA batteries & lanyard included.

#### CODE

MD1913

You can track your blood oxygen levels by inputting your readings manually into the Medicare LifeSense® App.



### EDAN Hand Held Pulse Oximeter

- Backlight control and automatic power-off function for power saving.
- SpO2 and Pulse rate measurement.
- Display numeric and waveform of SpO2 simultaneously. Up to 50 hours on 4 "AAA" batteries.
- Powerful data storage capacity.

#### CODE

EDA001



### EDAN SD1 Ultrasonic Foetal Pocket Doppler

Using conventional technology, EDAN SD1 breaks through the limitation of the probe angles as well as the depth of the foetal heart. The user can locate the foetal heart rate signal quicker and easier. The extremely low radiation level gives you a safer and more comfortable way during application.

#### CODE

EDA15

A photograph of a female doctor with grey hair, wearing a white lab coat and a stethoscope, checking the blood pressure of a male patient with dark skin. The patient is sitting on a black examination table, wearing a brown t-shirt and light-colored trousers. A blood pressure cuff is wrapped around his upper arm, and the doctor is holding the gauge. In the background, there is a red sharps container, a sink, and some medical equipment on a shelf. The image is partially covered by a large orange curved shape in the top right corner.

FURNITURE

## FURNITURE

Browse our range of premium quality healthcare furniture. Our range includes everything you need for your surgery and consultation room.



## Sunflower Gas Lift Stool

- Maximum user weight: 110kg (17st.).
- Anti-bacterial/anti-microbial vinyl.
- Luxury 65mm deep upholstery.
- 2 layer foam for support and comfort.
- Gas-assisted height adjustment.
- 42-58 height range.
- Five castor base for stability and manoeuvrability.
- Easy clean design and materials.

### CODE

SUN-STO2



## Sunflower Single Tier Couch Step - Rectangular

The Single Tier (Rectangular) Couch Step is lightweight and robust - an ideal addition for any surgery. They are stable and designed not to tip. The step has a non-slip surface and plastic feet to prevent movement.

- Maximum user weight 150kg (24st.)
- Anti-slip surfaces.
- Anti-slipping designs.
- Rubber feet prevent movement.
- Easy clean design and material.

### CODE

SUN-STP3



## Colenso Couch Black Upholstery

- 79cm (H) x 193cm (L) x 60cm (W) with paper roll holder.
- Maximum user weight: 165kg (26st.).
- Multi-position self locking backrest.
- 2 Layer/density foam for support & comfort.
- Luxury 65mm deep contoured upholstery.
- Anti-bacterial/anti-microbial vinyl.
- Paper roll holder supplied as standard.
- Durable silver epoxy coated frame.

### CODE

SUN-MEC2





### Sunflower Vista 50 Trolley

- Size: 101cm (H) x 56cm (L) x 52cm (W).
- Tilt and locking tray runner system.
- Trays integrate with all other Vista products.
- 100mm castors with front brakes.
- Optional tamper resistant locking door.
- Removable top shelf with integrated upstand.
- Easy clean design and materials.
- Supplied fully assembled.

#### CODE

SUN-MPT2



### Sunflower Vista 40 Trolley

- Tilt and locking tray runner system.
- Trays integrate with all other Vista products.
- Low level (80cm), ideal for under worktops.
- 100mm castors with front brakes.
- Optional tamper resistant locking door.
- Removable top shelf with integrated upstand.
- Stainless steel handle.
- Easy clean design and materials.
- Supplied fully assembled.

#### CODE

SUN-MPT41



### Sunflower Records Transfer Trolley

Designed to provide a complete solution for storing, distributing and transferring medical notes and files quickly to different locations. Records transfer trolleys are aesthetically pleasing, hygienically designed and easy to access, while being functional, strong, secure and easy to clean.

#### CODE

SUN-MRTTL1





## Sunflower Surgical Trolleys

- Mirror polished stainless steel shelves, trays and frame.
- Fully welded seamless frame.
- Fixed or removable/reversible shelves.

CODE	SIZE
SUN-STFW4	90cm H x 45cm L x 55cm W



- Anti-static shelf/tray protection buffers.
- Maximum individual tray load: 10kg.
- Supplied fully assembled.

CODE	SIZE
SUN-STFW7-RRFS2	86cm H x 75cm L x 45cm W



## Sunflower Four Panel Screen

- Provides a temporary or permanent privacy solution.
- Ideal for use as a hospital bed screen or as a general privacy screen.
- Hygienically designed using durable and easy to clean proven panel system.

CODE
SUN-MSF4



## Telescopic IV Pole

- 5 leg base tube with proxy grey finish.
- Sliding IV pole with chrome finish.
- Size: 2460mm x 1340mm.
- Floor space requirement: 580mm.
- Weight: 3kg.

CODE
ZHA01



SIMULATION

## SIMULATION

We offer a wide range of simulation products to facilitate effective training and development, such as our broad selection of training manikins, and demo dose blood and medication vials. Online training through the Fleming Medical Academy platform can be utilised for learning purposes.

Authorised distribution partner for

*Nasco*  
HEALTHCARE



## Life/form® Adult CRiSis™ Manikin

The Adult CRiSis™ Manikin comes complete with IV arm, blood pressure arm, compressible chest, standard airway management head, manual carotid pulse, constricted and dilated eye, shoulder IM injection, and hard carry case.

- Ideal for RCUK ALS, RCUK ILS, ACLS, Paramedic, EMT, and nursing training at every level.
- Measurements: 41" (H) x 19" (W) x 13" (D).
- Five year warranty.

### CODE

LF03890



## Simulaids WMD/CBRNE/DECON Full Body Trainer

The WMD/CBRNE/DECON Full-Body Trainer is a 5' 5", full-body trainer and combines standard intubation head and IV arm with a non-electronic body.

Basic life support services utilising advanced skills, like airway maintenance and IV administration, will appreciate this unit. The addition of a casualty simulation kit to this manikin broadens its use for MCIs.

Water won't damage this trainer, and you can initialise advanced life support treatment as you are flushing contaminants. (Please note: No ECG features).

- Size: 65" x 23" x 9".

### CODE

101-088FB



## Adult Choking Manikin

The Adult Choking Manikin with Carry Bag functions on the principle of clearing a blocked airway with forced air generated with the proper procedure, abdominal thrusts or back blows.

- Manikins include rib cage, xiphoid process, suprasternal notch to provide anatomical reference, foreign body for airway obstruction, shirt, and soft carry bag.
- Size: 30" x 17" x 10".

### CODE

101-1602



### Instructor's Economy Starter Package

This package is intended to give instructors the opportunity to replace older equipment or purchase start-up needs in the most economical way possible.

- 1 economy adult sani-manikin 4 pack with carry storage bag.
- 200 adult face shield lung systems.
- 1 sani-baby manikin 4 pack with carry storage bag.
- 200 baby face shield lung systems.
- Size: 25" x 23" x 19" .

#### CODE

100-2151



### EMT Casualty Simulation Kit

The EMT Kit has the largest variety of strap-on wounds for ease in getting the most impact with the least skill used during incident training preparation. Additional stick-on wounds and makeup accessories that allow for more detailed enhancement round out this kit.

- Size: 24" x 9" x 20" .

#### CODE

800-818



### Tension Pneumothorax Simulator

This model allows medical personnel to practice proper needle insertion which enables built-up air pressure to escape the thorax. The model accepts needle and thoracotomy in second intercostal space in the midclavicular line or fifth intercostal space in the midaxillary line.

- Complete with foot pump, replacement 2-piece pleural chamber, rib cage, and soft life-like skin.
- Soft carry bag included (procedural needle not included).
- Size: 19" x 19" x 9" .

#### CODE

150-1800



## Simulaids® Brad Adult™ CPR Manikin

The Simulaids® Brad Adult™ CPR Manikin features a realistic look and accurate anatomical landmarks, therefore, allowing for the proper application of CPR training techniques.

- Brad has a realistic head tilt with a chin lift to open the airway.
- Easily manipulated airway to simulate airway obstruction or choking situations.
- Long torso for even more realistic abdominal thrusts.
- Anatomical landmarks including sternum and rib cage plus substernal notch.
- Economically, individuals use CPR manikin mouth/nosepieces.
- User-friendly CPR manikin airway/lung system with a non-re-breathing valve that reduces cross-contamination.

### CODE

SI2801



## Simulaids® Brad™ VTA

The Brad™ VTA has the same look and feel of the life-like Brad™ CPR trainer, which has been a favourite of many instructors for years.

- Economical CPR manikin is made of soft, realistic vinyl plastic over polyurethane foam and rib structure to provide consistent high-quality training.
- Tested true for over a million compressions.
- Durable and easy to use.
- Includes carry bag, kneeling pads, 3 disposable lung airway systems, 3 channel design mouth/nosepieces, and a 9V battery.
- Latex-free.
- 1 year warranty.

### CODE

100-2865

### Features include:

- Long torso for realistic abdominal thrusts.
- Realistic head tilt and chin lift for opening airway.
- Anatomical landmarks including sternum, rib cage, and substernal notch.
- Airway manipulates easily to simulate airway obstruction or choking.

- VTA system combines observable chest rise, which verifies proper ventilation, with visual feedback rate, depth of compressions, and release on compressions.
- Rate light comes on with the correct rate.
- Second light comes on to show correct depth with proper release.
- Both lights on indicate the delivery of quality CPR.





### Life/form® GERi™ Complete Nursing Skills Manikin

GERi™ has an elderly appearance with skin wrinkles and folds. Complete life-like range of motion, realistic patient positioning, and non-pinching joints - moves like a person.

- Lightweight - approximately 28lbs.
- Full-size adult manikin - measures 58".
- Overall female appearance.
- Converts to male with removal of wig and attachment of male genitals.
- Ideal for OBRA required training.
- Five year warranty.
- Carry case not included.

#### CODE

LF04001U




### Life/form® KERi™ Complete Nursing Skills Manikin - Light

KERi™ has an ageless appearance. Complete life-like range of motion, realistic patient positioning, and non-pinching joints, allowing a life-like range at motion.

- Lightweight - approximately 28lbs.
- Full-size adult manikin - measures 58".
- Overall female appearance.
- Converts to male with removal of wig and attachment of male genitals.
- Ideal for OBRA required training.
- Carry case not included.

#### CODE

LF04020



## Life/form® TERi™ Androgynous Geriatric Patient Skills Trainer - Light

A comprehensive male/female elderly patient care trainer for medical simulation. Evolving from our GERi™/KERi™ full body manikins, it has a more realistic look and feel accurately representing the human anatomy. Weight is distributed to represent a real patient for lifting and carrying.

- Passive range of motion.
- Cardiovascular support.
- Intravenous, injection, subcutaneous injection sites.
- Nasogastric tube procedures.
- Enhanced aesthetics.
- Grooming and daily living assistance.
- Respiratory support.
- Indwelling male and female urinary catheterisation.
- Advanced skills.

### CODE

LF04301



## Life/form® TERi™ Androgynous Geriatric Patient Care Trainer - Light

A comprehensive male/female elderly patient care trainer for medical simulation. Evolving from our GERi™/KERi™ full body manikins, it has a more realistic look and feel accurately representing the human anatomy. Weight is distributed to represent a real patient for lifting and carrying.

- Passive range of motion.
- Enhanced aesthetics.
- Grooming and daily living assistance.
- Size: 5ft.
- Includes IV arm and BP arm.

### CODE

LF04300



### Obstetrical Manikin with Overlays

The Obstetrical Manikin is an anatomically correct pelvic model with full term newborn and placenta to give your trainees realistic practice in multiple techniques, and to learn the procedure for emergency childbirth.

- Disposable umbilical cords with clamps foetal.
- The life-like pelvic cavity has pronounced pelvic landmarks, located spinal column, angled birth canal, ileum, ischium, sacrum, sacro spinous ligaments, and greater sciatic notch.
- Easily replaceable extra vulva, powder to make simulated blood.
- A clear abdominal overlay to see positioning are all included.
- The unit is capable of demonstrating shoulder distocia, foetal attitude, foetal station, foetal lie, cord prolapse, compound presentation, and caesarean section incision locations can be discussed. The foetus provided is of a realistic newborn size with fontanelles and cranial sutures.
- A soft carry bag is included.
- Size: 20" x 13" x 16".

#### CODE

110-180



### Rescue Randy Combat Challenge 165lb Manikin

Developed for life-like patient handling, transportation, and extrication training. These manikins can be safely used in situations that are either too hazardous or too uncomfortable for human volunteers.

- Made of durable vinyl with 4,100 lbs test plastic-coated cables.
- Articulated joints and weight distribution according to human weight.
- Ideal for rescue and extrication training from pole top, confined spaces, collapsed buildings, smoke rooms, and ladder carry-down protocols.
- Turnout gear not included.
- Size: 55" x 27" x 13".

#### CODE

149-1435



## Life/form® Basic Child CRiSis™ Manikin

The Child CRiSis™ Manikin is a dramatic, state-of-the-art manikin designed for teaching individuals or groups life-saving techniques for children. Simulating a 5-year-old child, the Child CRiSis™ Manikin provides complete PALS training capabilities it combines all of the skill stations on one manikin, making training more realistic and giving students an appreciation for the smaller scale of the patient.

The Life/form® Child CRiSis™ Manikin is anatomically correct in detail. Realistic anatomy and landmarks include a mouth, tongue, oral pharynx, epiglottis, arytenoids, vocal cords, trachea, and oesophagus. In addition to the palpable and visual landmarks, the head, neck, and jaw are fully articulated, allowing a realistic jaw thrust feature.

- The length of the Child CRiSis™ Manikin is 46-3/4" and weighs 20lbs.
- Ideal for ACLS, nursing, paramedic and EMT training at every level.
- Hard carry case (LF03467U) sold separately.

### CODE

LF03765



## Life/form® Deluxe Child CRiSis™ Manikin with ECG Simulator

The Deluxe Child CRiSis™ Manikin is a dramatic, state-of-the-art manikin designed for teaching individuals or groups life-saving techniques for children. Simulating an 5 year old child, manikin provides complete PALS training capabilities it combines all of the skill stations on one manikin, making training more realistic and giving students an appreciation for the smaller scale of the patient.

- Blood pressure arm.
- CPR.
- Arrhythmia recognition.
- Femoral access.
- Intraosseous infusion.
- IV arm.

### CODE

LF03617



### STAT Baby Advanced

STAT Baby Advanced has all the features of a STAT Baby. It gives the instructor control of the physiological features of the baby to deliver challenging training to develop student assessment skills and easily tracks sessions for later review and debriefing.

Key functions include managing student information, easily create alternate scenarios or use those already existing, and efficient records management of the training, all managed through an easy iPad interface.

It satisfies every training need for post-neonatal care of paediatric patients, from introducing students to basic concepts to challenging medical staff in advance medical procedures.

- Size: 55" x 28" x 11".
- One year warranty.

#### CODE

101-8010



### Life/form® C.H.A.R.L.I.E Neonatal Resuscitation Simulator

**C:** Compressions, cardiac.

**H:** Heat compatible.

**A:** Airway, arterial access via the umbilicus.

**R:** Resuscitation.

**L:** Laryngeal mask capable for oxygen use or PPV.

**I:** Intravenous, intraosseous, intubation.

**E:** ECG.

C.H.A.R.L.I.E comes complete with:

- Airways, 2, standard and advanced.
- Baby powder.
- Blood.
- Carry bag.
- Defibrillation chest.
- Intraosseous bones and skin.
- IV bag.
- IV skin and veins for hand and foot.
- Male and female genitalia.
- Myelomeningocele.
- Needles.
- Omphalocele.
- Umbilicus.
- Instruction manual.
- Five year warranty.

#### CODE

LF01421





## Life/form® Advanced Venipuncture and Injection Arm - Light

The realism of the Life/form® Advanced Venipuncture and Injection Arm is truly amazing. The soft, flexible fingers are moulded separately with extreme attention paid to every detail - right down to the fingerprints. Flexion of the wrist helps students develop manipulation skills. The replaceable skin rolls as the veins are palpated, and a discernible pop is felt when entering the veins.

The Life/form® moulding process reproduces the fine details of the skin to make the arm look and feel alive. Valves in the veins can be seen and palpated on the skin surface. Simulated veins and skin are completely replaceable to keep this training arm looking and working like new.

- Complete replacement kits are available and easy to use.
- Five year warranty.

### CODE

LF01121



## Life/form® Adult Venipuncture and Injection Training Arm - Light

The Life/form® Venipuncture and Injection Training Arm is the ideal way to teach venipuncture and injection techniques including starting IVs. The skin rolls as you palpate the vein and the characteristic "pop" can be felt as the needle penetrates the vein.


Veins are accessible at the antecubital fossa, along the forearm, and at the back of the hand, making it possible to practice venipuncture at any of the common sites. Practice of intramuscular injection is equally realistic. A foam pad simulates muscle tissue, and the boundaries are identified by natural bone-like landmarks.

- Includes a Life/form® arm with skin, 22 gauge needle, butterfly set, infusion set, 3cc syringe, 12cc syringe, one pint Life/form® blood, fluid supply bags, instruction manual, and soft carry case. Fluid supply stand not included.
- Size: 35" x 12" x 6".
- Five year warranty.

### CODE

LF00698

**CAUTION:** Use only Life/form® replacements. **Do not** use conventional latex tubing. This will cause malfunction of the simulator.



## Life/form® Intramuscular Injection Simulator

A simulated bony structure is embedded in the torso and represents the superior end of the femur, or greater trochanter, the posterior superior and anterior superior iliac spines, and the sacrum. The bony structure provides palpable anatomical landmarks so students can identify proper injection sites. A section of the upper, outer quadrant of the left gluteal area is cut away to allow students to visualise the underlying structures. The gluteus medius and gluteus maximus muscles, sciatic nerve, and vascular structures are clearly shown.

Three types of intramuscular injections can be taught and practiced on this simulator: dorsogluteal, ventrogluteal, and vastus lateralis injections. As with all Life/form® replicas and simulators, the Intramuscular Injection Simulator is extremely realistic. Skin and muscle textures, as well as bone shape and position, closely resemble a live patient. As a result, perforation of the tissue with a needle duplicates the sensation of administering an actual injection.

- Comes complete with teaching guide, supply of syringes, and hard carrying case.
- 5 year warranty.
- Carton size: 22" x 16" x 10".

### CODE

LF00961



## Peritoneal Dialysis Simulator

For continuous ambulatory peritoneal dialysis. Designed to introduce patients, students, and nurses to the essentials of CAPD procedures and care. The realistic Life/form® torso offers a realistic method to demonstrate and practice peritoneal dialysis. Success with CAPD is dependent on the patient's following dialysis procedures with extreme care. This Life/form® simulator gives the patient an opportunity to gain the confidence necessary for prolonged success with CAPD.

- Complete with indwelling Tenckhoff catheter in a hard carry case.
- Dialysis materials not included.
- Five year warranty.

### CODE

LF01027



Demo Dose® Clear Ampule 1ml

Demo Dose® Clear Ampules are used to teach the proper handling and opening of medical ampules.

CODE	SIZE	PACK SIZE
PN01097U	1ml	100



Demo Dose® Clear Ampule 2ml

Teach students the proper method of ampule opening. Each glass ampule is physically scored at the neck for ease in opening.

CODE	SIZE	PACK SIZE
PN01095	2ml	10



Demo Dose® Clear Ampule 5ml

Demo Dose® Clear Ampules are used to teach the proper handling and opening of medical ampules.

CODE	SIZE	PACK SIZE
PN01098	5ml	100



Demo Dose® Simulated Red Blood Cells - O Negative

Nursing instructors can now teach this important procedure to students in the Skills Lab setting. Each bag features realistic barcoded labels and patient identification labels. Can be spiked, just like regular blood bags.

CODE	SIZE
PN01003	450ml

**WARNING:** For instrumental use only. Not for human or animal use.



### Demo Dose® Mini Vial

This versatile product is ideal for simulating preparation of intradermal and immunisations. Contains distilled water that has been tinted yellow.

CODE	SIZE
PN01248U	2ml



### Demo Dose® Nitr Spray

Realistically develop scenarios using nitroglycerin.

- Contains sterile water.
- Provides 400mcg simulated nitroglycerin.
- Contains 200 doses.

CODE	SIZE
PN01197	400mcg



### Demo Dose® Nitr Spray

Demo Dose® Lovenx is used for post op immobility scenarios to prevent DVT.

- Therapeutic class: anti-coagulant.
- Type: Pre-filled single-use syringe with safety needle.
- Dosage: 40mg/0.4ml.

CODE	SIZE
PN01094	0.4ml



### Demo Dose® Code Drug II Bundle Complete Set

Bundle includes key products to implement cardiac emergency scenarios. Each syringe contains distilled water.

CODE
PN01258

Set includes one of each medication listed below:

- Lidocain 2% (5ml).
- Atropin (10ml).
- Calcm Chlr 10% (10ml).
- Epinephrn (10ml).
- Sodium Bicarb 8.4% (50ml).
- Dextros 50% (50ml).
- Vasoprssn (1ml).
- Amiodron (3ml).
- Furosmide (10mg/ml).

**WARNING:** For instrumental use only. Not for human or animal use.



Demo Dose® Powder Vial

Now you can incorporate an antibiotic into your scenarios.

- Vial contains white powder for reconstitution for single use.

CODE	SIZE
PN01249	1g



Demo Dose® Powder Vial 10ml

Versatile product is ideal for simulating preparation of intradermal and immunisations. Contains distilled water that has been tinted yellow.

CODE	SIZE
PN01048U	10ml

**WARNING:** For instrumental use only. Not for human or animal use.



Demo Dose® Inject-Ed™ Injection Pad

Compact and economical injection task trainer that provides students with a realistic injection experience in the simulation lab. Only product currently on the market that allows you to inject fluid and easily drain out.

- Latex-free.
- IM, SQ, ID, and Z-track simulation.
- Air compatible.
- Resilient, durable, and self-healing surface is fluid (holds up to 20ml of fluid).
- Care instructions included.

CODE	SIZE
PN01195	3" x 1-1/2"



Demo Dose® Inject-Ed™ XL Injection Pad

Compact and economical injection task trainer that provides students with a realistic injection experience in the simulation lab.

- Latex-free.
- IM, SQ, ID, and Z-track simulation.
- Air compatible.
- Resilient, durable, and self-healing surface is fluid (holds up to 250ml of fluid).
- Care instructions included.

CODE	SIZE
PN01259	8" x 1-1/2" x 7"





FIRST AID

## FIRST AID

We are a leading supplier of First Aid Kits in Ireland and the UK, for healthcare, home, work and sport.

Our Medicare First Aid Kits are filled with premium quality products to ensure injuries are treated quickly and effectively.



Medicare First Aid Mini Kit  
1 - 3 PERSONS

- A more comprehensive kit that is great for storing in places like your car glovebox.

CODE	TYPE
MD6002	1 - 3 Persons



Medicare First Aid Travel Kit  
1 - 5 PERSONS

- Perfect for travelling.
- Water-resistant soft case with convenient belt loop and large easy-to-grab zip pulls.

CODE	TYPE
MD6003	1 - 5 Persons



Medicare First Aid Home Kit  
1 - 10 PERSONS

- The essential first aid home kit for life's little emergencies.

CODE	TYPE
MD6004	1 - 10 Persons



## Medicare First Aid Work Kit

### 1 - 10 PERSONS

The Medicare First Aid Work Kit 1 - 10 Persons is packed with all the essentials needed for any kind of emergency in the workplace.

- Perfect for the workplace.
- HSA approved work kit.

CODE	TYPE
MD6005	1 - 10 Persons



## Medicare First Aid Work Kit

### 10 - 25 PERSONS

Medicare First Aid Work Kit 10 - 25 Persons is packed with all the essentials needed for any kind of emergency in the workplace.

- Multi-compartment, durable hard case with ergonomic handle.
- Easy click on/off wall mounting bracket included.

CODE	TYPE
MD6006	10 - 25 Persons



## Medicare First Aid Work Kit

### 25 - 50 PERSONS

The Medicare First Aid Work Kit 25 - 50 Persons is packed with all the essentials needed for any kind of emergency in the workplace. Multi-compartment, durable hard case with ergonomic handle. Easy click on/off wall mounting bracket included, just grab and go. Includes essentials needed for everyday cuts and scrapes as well as emergency situations.

- Meets HSA Content Guidelines.

CODE	TYPE
MD6007	25 - 50 Persons



## Medicare Physiologix First Aid Sports Kit

The Medicare Physiologix First Aid Sports Kit is packed with all the essentials needed to treat a minor injury on the playing field. Suitable for all sports applications.

- Waterproof run-on bag.
- One main compartment with front and side pockets and water bottle holder.
- Contains a variety of dressings, tapes, bandages and cold therapy items for treating minor sports injuries.

CODE	TYPE
MD6008	Sports Kit



## Medicare CPR Face Shield

- Single use only.
- Designed to reduce the risk of cross-contamination.

CODE
MD1846



## Medicare Burn Gel

Medicare Sterile Burn Gel provides instant relief from burns, scalds and sunburn.

- Cools, soothes and moisturises the burn area and eases pain.
- Minimises trauma.
- Non-adherent - easily washed off with cool water or saline solution.

CODE	SIZE
MD77	3.5ml
MD78	50ml





### Medicare Burn Gel With Aloe Vera

Medicare Burn Gel is a water based gel with Aloe Vera for the first aid treatment of sunburns, minor burns, and scalds.

- Cooling effect on burns and pain relief.
- Transparent protective gel.
- Non-contact, easy and comfortable application.
- Can be used with the canister in any position.

CODE	SIZE
MD1220	100ml



### Medicare Burn Dressings

Medicare Burn Plasters are sterile hydrogel plasters suitable for the protection of minor and superficial burns.

CODE	SIZE
MD73	5 x 5cm
MD74	10 x 10cm
MD75	20 x 20cm
MD76	40 x 60cm

- Suitable for use on 1st, 2nd and 3rd degree burns and sunburns.
- Cools and soothes the burn and eases pain.

- Has the ability to absorb and dissipate heat.
- Promotes faster healing.
- Non-adherent and non-toxic.



### Medicare Burn Plasters 6's (Display of 16)

Medicare Burn Plasters are sterile hydrogel plasters suitable for the protection of minor and superficial burns.

- Provides comfort, cushioning and protection.
- Provides a bacterial barrier.
- Breathable.
- Contains Hydrogel which moisturises, cools and soothes, minimising trauma and helping to prevent infection.

CODE	PACK SIZE
MD190	6





## Medicare Spray Plaster

- Water resistant.
- Gradually disappears as the skin heals.
- Provides a transparent, breathable film.
- Ideal for awkward, hard-to-reach areas like the elbows, knuckles, toes or knees.
- Non-contact first aid that helps prevent infection.
- Ideal for use in situations where it is not recommended to touch the skin, for example schools.

CODE	SIZE
MD500	70ml



## Medicare Alcohol Prep Pads 50's

These soft, absorbent, non-woven alcohol wipes provide rapid-acting skin preparation in a sterile formula. Each pad is impregnated with 99.9% isopropyl alcohol. The non-woven swabs are manufactured to the highest bacteriological standards.

- Non-woven pads impregnated with 99.9% Isopropyl alcohol.
- Made from Biodegradable Crepe Paper.
- Sealed in a flexible foil sachet.
- Ideal for disinfection, antiseptic and topical microbial use.

CODE	PACK SIZE
MD1489	50



## Medicare Hot & Cold Reusable Pack

Recover faster with the natural healing power of hot and cold therapy.

- Heat Therapy: For effective relief from muscular aches and pains, joint stiffness, backache, strains and sprains.
- Cold Therapy: For fast relief from sprains, strains, bruises and swelling, headaches and minor injuries.
- Reusable.

CODE	PACK SIZE
MD602	10



## Medicare Sterile Saline Solution

Medicare Sterile Saline Solution is perfect for rinsing and cleansing away dust and dirt from wounds and eyes.

- Removes dust and dirt from wounds and eyes.
- Hypoallergenic.
- Sterile until the last drop.
- Can be used with the canister in any position.

CODE	SIZE
MD1221	100ml



## Medicare Emergency Eye Wash

Medicare Emergency Eyewash is a disposable eye irrigation solution. Emergency eye care is any type of eye care that is performed to stop immediate threat to vision.

- Disposable – for individual use.
- May be used daily as an eye wash or bath to cleanse the eyes.
- Use up to 6 times a day, as required.
- Supplied in boxes of 48 x 20ml vials.

CODE	SIZE	PACK SIZE
MD533	20ml	Box of 48



## Emergency Eye Wash

Emergency Sterile Eye Wash is used to minimise the effects of accidental exposure to chemicals. It is used for irrigation and cleansing of eyes.

- Can be used as much as needed.
- Contains sodium chloride solution.

CODE	SIZE
MD107	500ml



## Emergency Eye Wash Station

The Emergency Eye Wash Station provides on-the-spot decontamination and are a necessary back-up to minimise the effects of accident exposure to chemicals.

- Contains mirror and 2 x 500ml bottles of eye wash.
- Mirror is shatter-proof.

CODE	SIZE
CM11	2 x 500ml



Medicare Foil Rescue Blanket

- Medicare Foil Rescue Blanket prevents shock and is ideal for accident or emergency situations.
- Keeps in body heat and preserves body temperature.

CODE	SIZE
MD1845	Adult: 140 x 220cm



Medicare Triangular Calico Bandage

- Medicare Calico Triangular Bandage is a general purpose sling.
- Heard wearing, versatile and washable.

CODE	SIZE
MD1890	90 x 90 x 127cm



Medicare Skin Sanitiser Spray 70ml

Medicare Skin Sanitiser Spray is an alcohol-based skin sanitiser for protecting hands from germs and bacteria.

- Fast drying, non-sticky and rinse free.
- Anti-bacterial with 66% alcohol.
- Kills 99.9% of germs.

CODE	SIZE
MD587	70ml



Medicare Quick-Release Elasticated Tourniquet For Bleeding Control

- Quick Release Elasticated Tourniquet for the control of bleeding.
- Controls bleeding by stopping the flow of blood to the wound.
- One hand operation allows no pinching, jamming or catching.
- Quick release buckle provides instant removal.

CODE	SIZE
MD1840	100ml



Medicare Dressing  
Scissors Sharp/Blunt  
KT01



Medicare Dressing  
Scissors Sharp/Sharp  
KT02



Medicare Dressing  
Scissors Blunt/Blunt  
KT03



Medicare Lister  
Bandage Scissors  
KT04



Medicare First Aid  
Scissors  
KT05



Medicare Tough Cut  
Emergency Shears  
KT08

CODE	TYPE	SIZE
KT01	Dressing - Sharp/Blunt	14cm
KT02	Dressing - Sharp/Sharp	14cm
KT03	Dressing - Blunt/Blunt	14cm
KT04	Lister Bandage	14cm
KT05	First Aid	11cm
KT08	Tough Cut	18cm



## Medicare Fabric Plasters (100 Plasters)

- Hypoallergenic.
- Non-adherent wound pad.
- Highly breathable.
- 19 x 72mm Fabric Plasters.
- Value pack size of 100 plasters for regular use.

CODE	SIZE	PACK SIZE
MD207	19 x 72mm	100



## Medicare Blue Detectable Plasters (100 Plasters)

- 19 x 76mm Detectable Plasters.
- Made from PVC material.
- Easy visual detection in clean areas.
- Comes in an economic pack size of 100.

CODE	SIZE	PACK SIZE
MD209	19 x 76mm	100



## Medicare Waterproof Plasters (100 Plasters)

- 19 x 76mm Waterproof Plasters.
- Value pack size of 100 plasters for regular use.

CODE	SIZE	PACK SIZE
MD206	19 x 76mm	100



## Medicare Waterproof Spot Plasters (100 Plasters)

- Ideal shaped plaster for protecting the skin after injections.
- Hypoallergenic.
- Highly breathable.
- Non-adherent wound pad.
- Size: 2.2cm diameter.
- 100 pcs per pack.

CODE	SIZE	PACK SIZE
MD208	2.2cm	100



## **SECURE DRESSING FIXATION & COMPRESSION**

Our range of bandages offer firm directional support and protection for strains, sprains and minor burns.

# BANDAGES



## Mediform Conforming Bandage

The Mediform Conforming Bandage material is soft and highly adaptable. It can be used to hold wound dressings of all types in place with light compression effect.

- Made from polyamide and viscose.
- Comfortable, soft, lightweight, durable bandage.
- Conforms well to the body, ensuring patient comfort.

CODE	SIZE	PACK SIZE
MD0772	5cm x 4m	40
MD0773	8cm x 4m	40

CODE	SIZE	PACK SIZE
MD0774	10cm x 4m	20
MD0775	15cm x 4m	20



## BV Sterile Conforming Bandage

- Made from 41% cotton and 59% polyamide.
- Extremely conformable and effective at dressing retention.
- Individually wrapped.

CODE	SIZE	PACK SIZE
BA412-702	5cm x 4m	175
BA412-704	7cm x 4m	125
BA412-706	10cm x 4m	95
BA412-708	15cm x 4m	50



## Medicare Elastic Adhesive Bandage

Medicare Elastic Adhesive Bandages (EAB) offer firm directional support for strains, sprains and minor injuries. It is also used as a relieving and stretch compression bandage.

- Bidirectional easy tearing.
- Hypoallergenic skin friendly adhesive.
- Light and hard wearing.

CODE	SIZE	PACK SIZE
MIL2905	2.5cm x 4.5m	6
MIL2906	5cm x 4.5m	6

CODE	SIZE	PACK SIZE
MIL2907	7.5cm x 4.5m	6
MIL2908	10cm x 4.5m	6



CODE	SIZE	PACK SIZE
KO3482	5cm x 4.5m	12
KO3483	7.5cm x 4.5m	12

### Mediform Elastic Crepe Bandage

Medicrepe Elastic Crepe Bandage is a medium stretch bandage with constant elasticity and can be used for compression bandaging, securing dressings and devices, sprains and strains, protective bandaging, oedema control and partial immobilisation.

- Comfortable, soft, stretchable bandage, 100% cotton.
- Can be washed and sterilised.

CODE	SIZE	PACK SIZE
KO3484	10cm x 4.5m	12
KO3485	15cm x 4.5m	12



### BV Sterile Crepe Bandage

Sterile Crepe Bandage provides medium compression for light strains and sprains.

- Washable and reusable. 53% cotton, 38% polyamide and 9% polyurethane.
- Each bandage individually wrapped with clips.
- Double wrapped.

CODE	SIZE	PACK SIZE
BA424-703	10cm x 5m	65
BA424-705	15cm x 5m	50



CODE	SIZE	PACK SIZE
MD8976	6cm x 5m	24
MD8977	8cm x 5m	24

### Medicrepe Elastic Support Bandage

The Medicare Elastic Support Bandage comes with special clips to secure bandage once desirable level of compression has been achieved. Because of the high quality of fabric, there is superb elasticity.

- Cotton, polyamide and polyurethane high stretch bandage.
- Lightweight, fray resistant and durable.
- 50% cotton, 41% polyamide and 9% polyurethane.

CODE	SIZE	PACK SIZE
MD8978	10cm x 5m	12
MD8979	15cm x 5m	12

## BANDAGES



### BV Bandage Clip Fasteners

- Bandage Clip Fasteners with metallic ends for bandage fixation, can be used on all types of bandages.
- Avoids the use of adhesives.
- Can be easily repositioned.
- Allows for strong fixation.

CODE	COLOUR	PACK SIZE
BA470-002	Cream	30



### BV Elastic Cohesive Bandage

- Elastic Net Bandage is ideal for the fixation of dressings and other medical devices.
- Made from 70% polyamide and 30% polyurethane.
- Eliminates the use of clip fasteners or adhesives.

PACK SIZE	SIZE	PACK SIZE
BA4231-002	6cm x 4m	72
BA4231-003	8cm x 4m	72
BA4231-004	10cm x 4m	36



### PIC Self-Fix Elastic Cohesive Bandage

- Made of 85% cotton and 15% viscose.
- Zinc oxide in solvent base.
- Firm compression.
- Individually wrapped.
- Colour: White.

CODE	SIZE	PACK SIZE
AS30	10cm x 4m	6





### BV Compression Bandage

- High compression bandage with constant elasticity.
- Made from 85% cotton, 7% polyamide and 9% polyurethane.
- Washable, reusable and can be sterilised.
- Each bandage is individually wrapped with clips.

CODE	SIZE	PACK SIZE
BA416-018	8cm x 7m	180
BA416-019	10cm x 7m	140



### BV Tubular Gauze Bandage

- High compression bandage with constant elasticity.
- Made from 85% cotton, 7% polyamide and 9% polyurethane.
- Washable, reusable and can be sterilised.
- Each bandage is individually wrapped with clips.

CODE	SIZE	PACK SIZE
BA4212-001	6cm x 4m	486
BA4212-002	8cm x 4m	324
BA4212-004	12cm x 4m	216
BA4212-006	20cm x 4m	135



### BV Muslin Bandage

- White Open Weave Bandage is made from 100% cotton and 29 thread.
- For fixation of dressings and medical devices.
- Comes in a roll for a quick and easy dressing.
- Can be washed, sterilised and boiled for re-use.
- Individually wrapped.

CODE	SIZE	PACK SIZE
BA409-002	7cm x 5m	200
BA409-003	10cm x 5m	168





## BV Finger/Toe Tube Bandages

- Finger and Toe Tube Bandage made from 67% cotton and 33% viscose.
- Can be washed.
- Eliminates the use for clip fasteners or adhesives.
- Comes in a roll for a quick and easy finger dressing.
- clip fasteners or adhesives with this bandage.

CODE	SIZE	PACK SIZE
BA4211-100/1	Small	63
BA4211-101/1	Large	63



## BV Elasticated Viscose Stockinette

- Can be used as a dressing retention, without constriction or compression and to skin covering on any part of the body.
- Non-sterile, 2-way elastic bandage.
- Due to its manufacture techniques, the air can circulate between the threads, avoiding maceration.
- 90.5% viscose and 9.5% polyurethane.

CODE	SIZE	PACK SIZE
BA4215-132	7.5cm	12



## BV Stockinette

- Ribbed Stockinette is a 100% cotton, jersey, tubular ribbed bandage.
- Ideal for use as an underwrap for casts.
- Roll can be cut to the required length.

CODE	SIZE	PACK SIZE
BA4213-001	3cm x 25m	8
BA4213-002	5cm x 25m	20
BA4213-003	7cm x 25m	12
BA4213-004	10cm x 25m	8
BA4213-006	15cm x 25m	6



### BV Sterile Stockinette

- Ideal for use in orthopaedic surgery, traumatology and dermatology areas.
- 100% cotton, pre-cut and stitched.

CODE	SIZE	PACK SIZE
BA4216-001	5cm x 1m	6
BA4216-002	7.5cm x 1m	6
BA4216-003	10cm x 1.3m	6
BA4216-004	15cm x 1.4m	6
BA4216-005	20cm x 1.4m	6



### BV Elastic Net Bandage

- Elastic Net Bandage is ideal for the fixation of dressings and other medical devices.
- Made from 70% polyamide and 30% polyurethane.
- Eliminates the use of clip fasteners or adhesives.

CODE	SIZE	PACK SIZE
BA421-001	1.6cm x 25m	16
BA421-002	1.8cm x 25m	14
BA421-003	2.1cm x 25m	8
BA421-004	2.4cm x 25m	8
BA421-005	2.7cm x 25m	8
BA421-006	2.7cm x 25m	4
BA421-007	5.6cm x 25m	4



### BV Tubular Cotton Bandage

- Made of cotton.
- Protection of several parts of the body such as under orthopaedic or cast bandages.
- Sterilisable.
- Smooth and excellent dermophilic properties.
- Stretchable width-wise.

CODE	SIZE	PACK SIZE
BA4212-102	8cm x 25m	14
BA4212-105	15cm x 25m	8



## BV Elastic Tubular Support Bandage

- Elastic Tubular Support Bandages are made from 90% cotton, allowing for breathability.
- Excellent dermophil properties and very comfortable.
- Bandage can be doubled up.


CODE	SIZE	PACK SIZE
BA4214-001	4.2cm x 10m	8
BA4214-002	6.25m x 10m	6
BA4214-003	6.75cm x 10m	20
BA4214-004	7.5cm x 10m	20
BA4214-005	8.75cm x 10m	16
BA4214-006	10cm x 10m	20
BA4214-007	12cm x 10m	12
BA4214-009	21.5cm x 10m	10



## BV Elastic Tubular Bandage

- Made from 67% cotton and 33% viscose.
- Seamless sideways stretch-bandage is run-proof in any direction.

CODE	SIZE	PACK SIZE
BA4211-001	1.5cm x 20m	20
BA4211-002	2.5cm x 20m	12
BA4211-003	3.5cm x 20m	10
BA4211-004	5.5cm x 20m	6
BA4211-005	6.5cm x 20m	5
BA4211-006	7cm x 20m	5
BA4211-007	9cm x 20m	4
BA4211-008	12cm x 20m	3



### BV Tubular Bandage Applicators 6 Assorted Sizes

- Allows easy placing of elastic tubular bandages.
- Full set of applicators have six different diameters that cover all body applications from toes to thighs.
- Made from non-magnetic stainless steel.

CODE	SIZE	PACK SIZE
BA4210-001	See Below	1

#### Includes:

- (1) Size 1: 5 x 5 x 15cm
- (1) Size 2: 8 x 8 x 15cm
- (1) Size 3: 10 x 10 x 15cm
- (1) Size 4: 12 x 12 x 15cm
- (1) Size 5: 16 x 16 x 15cm
- (1) Size 6: 18 x 20 x 15cm



### BV Foam Collar and Cuff

- Soft foam covered sling/collar and cuff.
- Polyurethane foam.
- Covered with 100% cotton.
- 2 Pack.
- Can be cut to size.

CODE	SIZE	PACK SIZE
BA4218-001/1	5cm x 6m	12
BA4218-002/1	7.5 x 6m	8



### BV Gauze Rolls

BV Gauze Rolls provide fast-wicking action, superior aeration, excellent absorbency and reduces the risk of maceration.

- 20 thread/sq. cm.
- Folded 8ply.
- 100% cotton.

CODE	SIZE	PACK SIZE
BA402-001	80cm x 1m	246
BA402-002	80cm x 2m	120
BA402-003	80cm x 3m	56



## BV Sterile Orthopaedic Padding

- Orthopaedic padding is a non-woven sterile bandage of polyester layers, 85% polyester and 15% polyolefin.
- Ideal for use under a cast.

CODE	SIZE	PACK SIZE
BA420-702	7.5cm x 2.7m	240
BA420-703	10cm x 2.7m	168
BA420-704	15cm x 2.7m	108
BA420-708	10cm x 2.7m	168
BA420-709	15cm x 2.7m	108



## BV Non-Sterile Orthopaedic Padding

- Non-woven, non-sterile bandage of polyester layers.
- Ideal for use under a cast.
- Padding is sterilisable.

CODE	SIZE	PACK SIZE
BA420-001	5cm x 2.7m	288
BA420-002	7.5cm x 2.7m	192
BA420-003	10cm x 2.7m	216
BA420-004	15cm x 2.7m	96
BA420-005	20cm x 2.7m	104





ADVANCED  
WOUND CARE

## A SOLUTION FOR EVERY STAGE OF THE WOUND HEALING PROCESS

Our tried and trusted wound care products, which were designed alongside wound care industry experts, provide efficient, patient friendly remedies for conditions such as: acute wounds, chronic wounds, ulcers and burns.



### INJURY

### HAEMOSTASIS

### INFLAMMATORY

#### Stages of Wound Healing

##### Minutes


The wound is deep and bleeding heavily but no stitches are needed.



- Clean
- Blood coagulating into platelet plug

##### 0 - 3 Days

The wound is healing and debridement begins to occur. Exudate is starting to form.



- Surrounding skin may be red, hot, swollen
- White/yellow slough
- Minimal/moderate exudate

#### Treatment Objectives

- Cover and protect
- Absorb any remaining blood
- Moist wound healing

- Cover and protect
- Manage exudate
- Manage bacteria
- Aid removal of dead tissue
- Moist wound healing

#### Appropriate Dressing Category

- Primary Contact Layer
- Transparent Film
- Hydrocolloid
- Non-adherent Pad
- Surgical Dressing

- Foam
- Super Absorbent
- Alginate
- Hydrocolloid
- Amorphous Hydrogel
- Antimicrobial if at high risk of local infection

#### Advanced Woundcare Dressings



Including traumatic wounds, pressure ulcers, venous leg ulcers, diabetic foot ulcers (under medical/nursing supervision)

PROLIFERATION	EPITHELISATION	MATURATION
<b>3 - 24 Days</b> This is the reconstruction/repair phase. The wound fills from bottom up with new granulation tissue.	<b>3 - 24 Days</b> Still in the reconstruction/repair phase. New skin is forming over the wound and is very delicate.	<b>24 Days - 2 Years</b> The wound has healed but the skin is still delicate and therefore is a risk of reopening of the scar.
		
<ul style="list-style-type: none"> <li>• Red granulation tissue in wound bed</li> <li>• As wound fills with granulation tissue, pink/white epithelial tissue migrates from wound edge</li> <li>• Minimal or no exudate</li> </ul>	<ul style="list-style-type: none"> <li>• Pink/white epithelial tissue covers the wound</li> <li>• Very delicate</li> <li>• No exudate</li> </ul>	<ul style="list-style-type: none"> <li>• Wound is closed but not fully healed</li> </ul>
<ul style="list-style-type: none"> <li>• Cover and protect new tissue</li> <li>• Absorb exudate</li> <li>• Moist wound healing</li> </ul>	<ul style="list-style-type: none"> <li>• Cover and protect new epithelial tissue</li> <li>• Moist wound healing</li> </ul>	<ul style="list-style-type: none"> <li>• Cover and protect if at risk of new damage</li> </ul>
<ul style="list-style-type: none"> <li>• Foam</li> <li>• Hydrocolloid</li> </ul>	<ul style="list-style-type: none"> <li>• Transparent Film</li> <li>• Hydrocolloid</li> <li>• Non-adherent Pad</li> <li>• Surgical Dressing</li> <li>• Foam</li> </ul>	<ul style="list-style-type: none"> <li>• Transparent Film</li> <li>• Hydrocolloid</li> <li>• Non-adherent Pad</li> <li>• Surgical Dressing</li> </ul>
 	   	 



## Mediporex Adhesive Surgical Dressing

- Adhesive.
- Sterile and individually packed.

CODE	SIZE	PACK SIZE
MD1150	6 x 7cm	50
MD1151	9 x 10cm	50
MD1152	9 x 15cm	30

CODE	SIZE	PACK SIZE
MD1153	9 x 20cm	30
MD1154	9 x 30cm	30
MD1155	7 x 8cm	50



## Mediporex Waterproof Surgical Dressing

- Sterile.
- Coated with a vapour permeable film which provides a showerproof barrier.

CODE	SIZE	PACK SIZE
MD1159	7 x 8cm	50
MD1160	5 x 7cm	50
MD1161	10 x 10cm	50

CODE	SIZE	PACK SIZE
MD1162	10 x 15cm	30
MD1163	10 x 20cm	30



## Mediflex Silicone Foam Dressing

- Waterproof and bacteria proof backing.
- Sterilised and individually packaged.

CODE	SIZE	PACK SIZE
MD1164	6 x 7cm	5
MD1166	15 x 5cm	5



## Mediflex Bordered Silicone Foam Dressing

- Silicone adhesive border surrounding wound pad.
- Waterproof and bacteria proof backing.
- Sterilised and individually packaged.

CODE	SIZE	PACK SIZE
MD1171	7.5 x 7.5cm	5
MD1170	10 x 10cm	5



## Medifoam Heel Foam Dressing

- Specifically developed to address wounds on the heel or similar anatomical areas such as the heel or elbow where a standard dressing is not suitable.
- Non-adhesive.

CODE	SIZE	PACK SIZE
MD1384/1	12 x 13cm	3

## Medifoam Adhesive Foam Dressing

- Self-adhesive border dressing.
- Waterproof and bacteria proof backing.
- Sterilised and individually packaged.

CODE	SIZE	PACK SIZE
MD1388/1	10 x 10cm	5
MD1389/1	15 x 15cm	3

## Medifoam Non-Adhesive Foam Dressing

- Non-adhesive.
- Waterproof and bacteria proof backing.
- Sterilised and individually packaged.

CODE	SIZE	PACK SIZE
MD1386/1	10 x 10cm	5





## Medigauze Gauze Swabs

- Medigauze swabs are suitable for controlling bleeding.
- Highly absorbent material, 8 ply 100% surgical grade absorbent cotton gauze.
- Low fibre release reducing the risk of wound contamination.
- Soft material that is gentle on the skin.

CODE	SIZE	PACK SIZE
MD303	5 x 5cm	100
MD304	7.5 x 7.5cm	100
MD305	10 x 10cm	100



## Mediswab Non-Sterile Non-Woven Swabs

- Non-Sterile Swabs are used for many medical purposes.
- Suitable for use with all antiseptic solutions.

CODE	SIZE	PACK SIZE
MD1304/1	5 x 5cm	20
MD1303	7.5 x 7.5cm	20
MD1302	10 x 10cm	20



## Mediswab Non-Sterile Non-Woven Swabs

- Non-Sterile Swabs are used for many medical purposes.
- Suitable for use with all antiseptic solutions.

CODE	SIZE	PACK SIZE
MD300	5 x 5cm	100
MD301	7.5 x 7.5cm	100
MD302	10 x 10cm	100



## Medipad Non-Adherent Wound Pad

- Sterile.
- Suitable for use with all antiseptic solutions.

CODE	SIZE	PACK SIZE
MD1350/1	5 x 5cm	10
MD1351/1	10 x 10cm	10



## Medisite Transparent Island Dressing

- Highly absorbent: non-adherent central wound pad.

CODE	SIZE	PACK SIZE
MD1530	6 x 7cm	20
MD1531	10 x 12cm	20



## Medisite Transparent Thin Film Dressing

- Transparent film allows healing process to be observed without removing the dressing.

CODE	SIZE	PACK SIZE
MD1543	15 x 26cm	20



## Medisite Transparent Dressing with Grid

A sterile, transparent, polyurethane, bacteria-proof and waterproof film, coated with a hypoallergenic adhesive and supplied with a wound mapping grid carrier.

- Wound mapping grid allows simple but accurate measurements of the healing process to be made.

CODE	SIZE	PACK SIZE
MD1540	6 x 7cm	20



## Medistrip Wound Closure Strips

- Sterile.
- Suitable for use with all antiseptic solutions.

CODE	SIZE	PACK SIZE
MD1400/1	6 x 76mm	25
MD1401/1	6 x 102mm	25
MD1402/1	10 x 102mm	25



## Meditull Paraffin Gauze Dressing

- The dressing is used as a primary wound contact layer.

CODE	SIZE	PACK SIZE
MD1355	5 x 5cm	10
MD1354/1	10 x 10cm	10



## Medigel Super Absorbent Dressing

- Locks fluid cleanly away from the wound and skin.
- Fluid repellent backing layer minimises fluid strike-through.
- Reduction of dressing changes leading to increased comfort for the patient.

CODE	SIZE	PACK SIZE
MD1408	15 x 20cm	10



## Medifixx Dressing Tape

- Breathable. The porous structure of the non-woven fabric ensures air and vapour permeability.
- Skin-friendly, water-based adhesive.
- Designed to be gentle to the skin, but at the same time provide secure fixation.
- Comfortable, stretchable fabric that conforms well to body contours.

CODE	SIZE
MD1111	5cm x 10m
MD1112	10cm x 10m
MD1113	15cm x 10m



## Medipad Absorbent Pads

- Non-sterile.
- Suitable for use on moderate or heavy exuding wounds.
- Pad is made up of 4 layers.
- Extra absorbent.
- Breathable.
- Waterproof backing.

CODE	SIZE	PACK SIZE
MD142	20 x 20cm	25



## BV Ribbon Gauze

- 100% cotton., 22 thread/sq. cm.

CODE	TYPE	SIZE	PACK SIZE
BA425-001	Non-Sterile	1.25cm x 5m	24
BA425-002	Non-Sterile	2.5cm x 5m	12
BA425-016	Non-Sterile	1.25cm x 10m	24
BA425-017	Non-Sterile	2.5cm x 10m	12
BA425-502/1	Sterile	2.5cm x 5m	30
BA425-503/1	Sterile	5cm x 5m	26



## BV Sterile Absorbent Pads

- Sterile Absorbent Pads are suitable for use on moderate-heavy exuding wounds.
- Individually sealed and sterilised by steam or gas.
- Pad is made up of 4 layers.
- Waterproof backing.


CODE	SIZE	PACK SIZE
BA485-501	10 x 10cm	60
BA485-502	10 x 15cm	30
BA485-503	10 x 20cm	30
BA485-504	10 x 30cm	30
BA485-507	20 x 20cm	30
BA485-510	20 x 40cm	25



## BV Non-Sterile Absorbent Pads

- Suitable for use on moderate or heavy exuding wounds.
- Pad is made up of 4 layers.
- Waterproof backing.

CODE	SIZE	PACK SIZE
BA485-010	20 x 40cm	25



## BV Sterile Non-Adherent Dressing

- Made up of a non-adherent layer, an absorbent fleece, a fusible powder and a water repellent fabric layer.
- Absorbs up to 20 times its own weight in exudate.
- Individually wrapped.
- Steam sterilised at 121°C.

CODE	SIZE	PACK SIZE
BA487-001/1	5 x 5cm	150
BA487-002/1	7.5 x 7.5cm	150
BA487-003/1	10 x 10cm	150

CODE	SIZE	PACK SIZE
BA487-004/1	10 x 32cm	150
BA487-503/1	30 x 40cm	40



## BV Eye Pads Sterile

- Post-operative eye care.
- For the treatment of discharging or dry wounds.
- Excellent padding effect.
- Very soft and absorbent.
- 100% cotton.
- Approx. weight: 400g/sq. m.
- Individual peel-open pack.

CODE	SIZE	PACK SIZE
BA4372-001	5.6 x 7cm	60



## Medisilk Silk Tape

Medisilk Silk Tape is ideal for the fixation of dressing pads and medical devices. This is a silk finished, non-stretch, water repellent, tape with a strong, hypoallergenic adhesive. The porous, hypoallergenic tape makes it ideal for use on patients with sensitive skin, as it allows the skin to breathe through it.



CODE	SIZE	PACK SIZE
MD752	1.25cm x 5m	18
MD753	2.5cm x 5m	12
MD766	5cm x 5m	6
MD753/1	2.5cm x 9m	12

MD753/1: For Clinical use only.

9m Tape does not come with a plastic cover.

## Medilite Paper Tape

Medilite Paper Tape is a gentle, breathable paper tape. It is ideal for fixation of dressings where frequent taping/untaping is required and is perfect for use on patients with sensitive, delicate skin.



CODE	SIZE	PACK SIZE
MD754	1.25cm x 5m	18
MD755	2.5cm x 5m	12
MD767	5cm x 5m	6
MD755/1	2.5cm x 9m	12

9m Tape does not come with a plastic cover.

## Meditrans Transparent Tape

Transparent, perforated polyethylene film with bidirectional tear. Easily tears into very thin strips for securing tubing and dressings that need to be monitored.

- Ventilated - Allows skin to breathe.
- Skin friendly, due to hypoallergenic adhesive.



CODE	SIZE	PACK SIZE
MD758	1.25cm x 5m	18
MD757	2.5cm x 5m	12
MD757/1	2.5cm x 9m	12



## Meditape Zinc Oxide Stretch Tape

Highly conformable Zinc Oxide Tape that stretches in all directions. For compression applications and securing dressings in challenging areas. A stronger, more adhesive tape that stays fixed for longer. Ideal for strapping and taping in first aid and sports applications.

CODE	SIZE	PACK SIZE
MD756	2.5cm x 4m	12



## Meditape Zinc Oxide Tape

Meditape Zinc Oxide Tape is ideal for strong fixation of non-occlusive dressing material and medical devices. A stronger, more adhesive tape that stays fixed for longer. A micro-perforated, taffeta tape.

- Bidirectional easy tearing.
- Hypoallergenic skin friendly adhesive.
- Beige colour.

CODE	SIZE	PACK SIZE
MD750	1.25cm x 5m	18
MD751	2.5cm x 5m	12
MD765	5cm x 5m	6

## PROCEDURE PACKS

We offer a range of sterile medical sets for the treatment of wounds, and minor and emergency medical procedures. All the elements included in the procedure packs have been carefully selected and have been developed alongside our customers in response to specific procedural requirements.



## Mediset Dressing Set No.1

ITEM	QTY
Wrapping Field: 40 x 40cm	1
Tray: 15 x 7.5cm	1
Gauze Swabs: 4 x 5cm	2
Gauze Balls: Size D	3
Forceps	1

### CODE

MD400



## Mediset Dressing Set No.2

ITEM	QTY
Wrapping Field: 50 x 50cm	1
Tray: 18 x 9.5cm	1
Gauze Swabs: 4 x 5cm	4
Gauze Swabs: 7.5 x 7.5cm	2
Gauze Balls: Size D	2
Forceps	2

### CODE

MD401



## Mediset Dressing Set No.4

ITEM	QTY
Wrapping Field: 60 x 60cm	1
Tray: 18 x 9.5cm	1
8ply Compresses: 10 x 10cm	2
Gallipot	1
Gauze Balls 5.0cm	2
Green Tweezer	1
Waste Bag	1
Pair of Gloves	1

### CODE

MD403



## Mediset Catheterisation Set No.1

### CODE

MD405

ITEM	QTY
Green Tweezer	1
Kidney Tray	1
Wrapping Drape: 60 x 60cm	1
Gallipot: 250ml	1

ITEM	QTY
Waste Bag	1
Gauze Balls 5cm	4
1 Pair of Gloves	1



## Mediset Catheterisation Set No.3

### CODE

MD407

ITEM	QTY
Forceps: Blue/Green	2
Kidney Tray	1
Wrapping Drape: 60 x 60cm	1
Cut Drape: 60 x 60cm	2
Folded Drape: 60 x 60cm	1
Gallipot: 250ml	1

ITEM	QTY
Waste Bag	1
Gauze Balls 5cm	2
Gauze Balls 3cm	5
Tray: 250ml	1
8ply Compress: 10 x 10cm	2
Pair of Gloves	1



## Mediset Sterile Mouthcare Set

### ITEM

### QTY

Cleaning Sticks: PS Stick with Foam Head	1
Tray: 18 x 9.5cm x 2.5cm	5
Plastic Tongue Depressor	1
Wrapping Field: 50 x 50cm	1

### CODE

MD411





## BV Sterile Suture Pack No.2

ITEM	QTY
Wrapping Field: 50 x 50cm	1
Tray: 18 x 9.5cm	1
Gauze Swabs: 7.5 x 9cm	5
Gauze Balls: Size D	5
Forceps	1
Suture Blade	1
Pair of Gloves	1
Waste Bag	1
Clamp	1

## CODE

BA490-021



## BV Emergency Delivery Set

The kit is presented in carton box (27 x 20 x 11 cm) that has 2 bags, marked A for routines before delivery and B for routines after delivery. Each box also includes 8 disinfection wipes with alcohol and 2 pairs of surgical gloves.

## CODE

BA4901-060

## SET A - BEFORE DELIVERY

ITEM	QTY
Non-Woven Gown	1
Polyethylene Apron	1
Protection Sheet: 120 x 150 cm	1
Non-Woven Swabs: 10 x 10cm	5
Hand Towels: 35 x 36cm	4
Face Masks with Elastics	2
Waste Bag	1

### SET B - AFTER DELIVERY

ITEM	QTY
Baby Blanket: 90 x 80cm	1
Suction Bulb	1
Maternity Pad: 30 x 10cm	1
Non-Woven Swabs: 10 x 10cm	6
Waste Bag	1
Umbilical Cord Clamp	2
Hand Towels: 35 x 36cm	4
Laparotomy Swab: 30 x 45cm	1
Umbilical Cord Scissors	1



SURGICAL

## SURGICAL

Discover our surgical product range which includes dressing sets to facilitate the effective management of acute wounds such as lacerations, incisions, abrasions and invasive lines. We also offer a range of theatre gowns and scrub suits.



## BV Theatre Gowns Sterile Large Blue

- Non-woven Sterile Theatre Gown with long sleeves.
- Attached ties on back and neck.
- Knitted cuffs reinforced in chest and arms.
- Inside and outside belt sown to gown.
- Comes with 2 towels.
- Colour: Blue.
- Size: Large.
- Latex-free.

CODE	SIZE	PACK SIZE
BA462-602/1	Large	1



## Premier Protect Barrier Gown Standard

This gown is manufactured from soft durable polyethylene material of 1 piece construction. Impervious to fluids, providing an extremely effective barrier against cross infection.

- Tear away halter neck allows easy removal of the gown.
- Independently tested as 100% effective against strike through by staphylococcus aureus (MRSA).
- Thumb loops prevent leakage between glove and sleeve.

CODE	COLOUR	PACK SIZE
PM5530	Blue	20
PM5550	Green	20



## BV Non Sterile Scrub Suit

- Polypropylene non-woven.
- Approx. Weight: 38g/sq.m.
- Tunic with V collar, short sleeve and 1 pocket. Pants with elastic on the waist.
- Individual plastic bag (1 tunic and 1 pant) and transport carton.

CODE	SIZE	PACK SIZE
BA462-715	Large	1 Tunic, 1 Pant
BA462-716	Extra Large	1 Tunic, 1 Pant



### BV Examination Gown

- Polypropylene non-woven.
- Approx. weight: 40g/sq.m.
- Trespass model with short sleeve. Closes with 6 velcros in front. With two possibilities to close (tight and loose). Open in armpits.
- One size, fits all.
- Individual plastic bag and transport carton.
- Colour: Green.

CODE	SIZE	PACK SIZE
BA462-004	96 x 136cm	1 Gown



### BV Impervious Visitor Gown

- Laminated of non-woven polypropylene with polyethylene.
- Open in back, knitted cuffs, attached ties on neck and on waist.
- One size, fits all.
- Plastic bag and transport carton.

CODE	SIZE	PACK SIZE
BA462-009	OSFA	10



### BV Sterile Surgical Adhesive Drapes

- Sterile Surgical Adhesive Drapes consist of a laminate of green non-woven 100% Viscose and Polypropylene Film.
- Individual Peel Open Packs.

CODE	SIZE	PACK SIZE
BA488-203/1	75 x 90cm	25 per Pack
BA488-103/1	75 x 90cm	30 per Pack
BA488-101	45 x 75cm	50 per Pack



### Disposable PE Overshoes

- Disposable PE overshoes.
- Multi-purpose shoe covers.
- The elastic keeps the cover securely in place.

CODE	COLOUR	PACK SIZE
W025	Blue	100





## BV X-Ray Sterile Abdominal Swabs 12ply with Tape

- Sterile 12ply X-Ray Detectable Laparotomy Sponges (with tape) are made from 13 thread bleached 100% pure cotton.
- X-ray thread is monofilament type, with at least 60% of barium sulphate, one thread by layer.
- Machine folded with superior finish.
- Sewing thread also 100% cotton and unbleached without any optical agents.
- 100 pcs (20 x 5 packs).

CODE	SIZE	PACK SIZE
BA492-059	30 x 30cm	100



## BV X-Ray Detectable Sterile Swabs 32ply

- Made from 13 thread bleached 100% pure cotton.
- Sewing thread is 100% cotton unbleached without optical agents.
- X-Ray thread is monofilament with at least 60% barium sulphate.
- Folded with tucked in edges to make sure no loose fibres can fall into wound bed.

CODE	SIZE	PACK SIZE
BA492-202	10 x 10cm	200
BA492-209	5 x 5cm	550



## BV X-Ray Non-Sterile Gauze Swabs

- Non-sterile 12ply X-ray detectable gauze swabs.
- Sewing thread is 100% cotton, unbleached without optical agents.
- X-ray thread is monofilament with at least 60% of barium sulphate.
- Individually wrapped.

CODE	SIZE	PACK SIZE
BA428-053	5 x 5cm	100





### BV Tracheotomy Sponge

- Sterile Tracheotomy Sponges consist of a laminate of 3 layers: Non-adherent fleece, absorbent fleece and water repellent viscose.
- Expandable hole to accommodate tubes of diameter 12-16mm.
- Packed in a sterile peel-open pack.

CODE	SIZE	PACK SIZE
BA4371-002	8 x 9cm	150



### BV X-Ray Sterile Detectable Laparotomy Sponges 12ply

- 12ply sponges (without tape).
- Made from 13 thread bleached 100% pure cotton.
- X-ray thread is monofilament type, with at least 60% of barium sulphate, one thread by layer.
- 100 pcs (20 x 5 packs).

CODE	SIZE	PACK SIZE
BA492-058	30 x 30cm	100



### BV Kidney Trays

- Sterile and graduated.
- 300, 800 and 1500ml kidney tray with a pour spout for easier transference of liquids.
- Individual peel-open pack and transport carton.

CODE	SIZE
BA489-510	700ml
BA489-516	300ml
BA489-520	1500ml



## BV X-Ray Detectable Sterile Gauze Balls

- Made from 20 thread/sq. cm.
- PP X-ray thread.
- Packed in sterile peel-open packs.
- 125 pcs (25 x 5 packs).

CODE	SIZE	PACK SIZE
BA475-039/1	Size F: 50mm	125



## BV Non-Sterile Gauze Balls

- Non-sterile gauze balls made from 17 thread bleached 100% pure cotton.
- All edges turned in to avoid loose threads.
- Elastic ring can stand steam sterilisation up to 135°C.

CODE	SIZE	PACK SIZE
BA429-005	Size E: 30 x 40cm	100



## BV Sterile Gauze Balls without Elastic Ring

- Sterile gauze balls without elastic ring is made from 20 thread/sq. cm.
- 100% cotton.
- Edges folded inside to avoid loose threads.
- Packed in sterile peel-open packs.
- 100 pcs (20 x 5 packs).

CODE	SIZE	PACK SIZE
BA475-030	Size H: 55 x 70cm	100



## BV Non-Sterile Non-Woven Ball with Elastic Ring

- Made from 70% viscose and 30% polyester.
- Latex-free.
- Edges turned inside to avoid loose threads.
- Elastic ring can stand steam sterilisation up to 135°C.


CODE	TYPE	SIZE	PACK SIZE
BA440-003	Non-Sterile	Size F: 17 x 20cm	100
BA483-500	Sterile	Size E: 30 x 35cm	800



## INFUSION

## INFUSION

Our infusion range includes a variety of needles and syringes that enable the safe passing of fluids into the bloodstream and general phlebotomy use.



## PIC Hypodermic Needle 100's

- Hypodermic Needles are non-pyrogenic and non-toxic disposable sterile needles.
- Wide range of gauges allows its uses in a variety of treatments including: interdermal, subcutaneous, intramuscular, intravenous and intraarticular.
- Each needle has a luer hub and is individually packed in a peel-open blister pack.

CODE	COLOUR	SIZE	PACK SIZE
AN02	Pink	18g x 1.5"	100
AN03	Cream	19g x 1.5"	100
AN04	Yellow	20g x 1"	100

CODE	COLOUR	SIZE	PACK SIZE
AN05	Green	21g x 1.5"	100
AN08	Yellow	30g x 0.5"	100



## PIC Hypodermic Needle 100's

CODE	SIZE	COLOUR	PACK SIZE
AN36	25g x 0.625"	Orange	100
AN38	23g x 0.625"	Orange	100

- Hypodermic Needles are non-pyrogenic and non-toxic disposable sterile needles.
- Wide range of gauges allows its uses in a variety of treatments including: interdermal, subcutaneous, intramuscular, intravenous and intraarticular.

- Each needle has a luer hub and is individually packed in a peel-open blister pack.
- Sterilised with ethylene oxide, have thin walls and high penetration.



## PIC Hypodermic Syringe 1ml - Insulin Syringe 160's

- Size: 1ml.
- 160 pcs per box.

CODE	SIZE	PACK SIZE
AN40	1ml	160



## Boscarol Solution Administration Set

- Latex-free.
- Single use only.

CODE	TYPE
BOS42	1 Port

## i Syringe Tip Types

### Luer Slip (or Slip Tip)

Allows a needle to be quickly and conveniently pushed straight on to the end of the tip (if a needle is required). So if you are looking for speed of connection between the needle and syringe, the Luer Slip will suit you best.

### Luer Lock

Allows a Luer Lock needle to be screwed onto the tip and securely locked in place ensuring no liquids can leak. So if you are looking for a more secure fitting with less chance of leakage - and therefore more accurate measurement and application - choose the Luer Lock tip.

### Luer Slip Eccentric

Used when you need to inject something parallel to the skin. The off-centre tip is in line with the outer side of the barrel, so is perfect for injecting into a vein or artery.

### Catheter

Commonly used when injecting via tubing, flushing wounds or when using a slip tip needle that is too large to use standard slip tip syringe.



- Ideal for drawing of biological matter for diagnostic purposes (according to the selected needle).
- Disposable and sterile.
- Latex Free.

## PIC Luer Slip Syringe without Needle

CODE	SIZE	PACK SIZE
AN20	2.5ml luer slip	100
AN11	5ml luer slip	100
AN41	10ml luer slip	100

- Made from high transparency polypropylene.
- Triple contact rubber plunger.
- Packed in individual pel-open blister packs with safety device.



## PIC Luer Slip Eccentric Syringe without Needle

- Ideal for administration/infusion of solutions.
- Made from high transparency polypropylene.
- Triple contact rubber plunger.
- Disposable and sterile.
- Each pack of syringes has a magnifying glass to help measure dosage.

CODE	SIZE	PACK SIZE
AN12/3	10ml luer slip eccentric	100
AN15	50ml luer slip eccentric	50
AN35	20ml luer slip eccentric	100





## PIC Luer Lock Syringe without Needle

- Suitable for administration/infusion of solutions.
- Ideal for drawing of biological matter for diagnostic purposes (according to the selected needle).
- Latex-free.
- Disposable and sterile.
- Made from high transparency polypropylene.
- Triple contact rubber gasket and plunger to prevent accident leakage during aspiration.

CODE	SIZE	PACK SIZE
AN17	50ml luer lock	50
AN30	20ml luer lock	100
AN37	10ml luer lock	100
AN20	2.5ml luer lock	100



## PIC Catheter Tip Syringe without Needle

- Catheter Tip Syringe (without needle) is suitable for administration/infusion of solutions.
- Ideal for drawing of biological matter for diagnostic purposes (according to the selected needle).
- Disposable and sterile.
- Made from high transparency polypropylene.
- Latex-free.
- Packed in individual peel-open blister packs with safety advice.

CODE	SIZE	PACK SIZE
AN39	50ml catheter	50



## Pic 2-Way IV Cannula Teflon Catheter

- Suitable for peripheral intravenous procedures for therapeutic and diagnostic purposes.
- Sterile, disposable, non-pyrogenic and non-toxic.
- Transparent chamber has hydrophobic filter for immediate visualisation of blood back-flow.
- Perfect needle alignment between needle point and catheter conical sections.
- Colour-coded needles.

CODE	SIZE	COLOUR	PACK SIZE
AN24	18g x 45mm	Green	50
AN43	22g x 25mm	Blue	100



## OXYGEN MASKS & AIRWAYS

## OXYGEN MASKS & AIRWAYS

Discover our high-quality range of oxygen masks and airways equipment designed for effective patient treatment and care.



### Wellead Oxygen Mask with Tubing

- Adult oxygen mask with tubing is made from clear soft PCV.
- Adjustable nose clip assures comfortable fit.
- Latex-free.
- 210cm oxygen supply tubing can be configured (medium concentration).
- Batch marking for traceability (applicable 40-50% in clinical use regarding oxygen concentration of masks).

CODE	TYPE
WL331003	Adult



### Boscarol Oxygen Mask with Tubing

- Adult oxygen mask with tubing is made from clear soft PCV.
- Adjustable nose clip assures comfortable fit.
- Latex-free.
- 210cm oxygen supply tubing can be configured (medium concentration).
- Batch marking for traceability (applicable 40-50% in clinical use regarding oxygen concentration of masks).

CODE	TYPE
BOS81	Adult



### Wellead Disposable Adult Nebuliser Kit with 6ml Jar

- Nebuliser with aerosol mask is made from clear, soft and non-toxic PVC.
- Helix connector.
- Latex-free.
- Batch marking for traceability.
- CE class IIA.
- Includes nebuliser mask (6mm), jar and tubing.
- 7ft tubing.

CODE	TYPE
WL343002	Adult
WL343004	Paediatric



## Wellead Non-Rebreather Mask

- Made from non-toxic PVC.
- Latex-free.
- Bench marking for traceability.
- CE class IIA.

CODE	TYPE
WL333002	Adult
WL333004	Paediatric



## Wellead Yankauer Handle without Vent

- Made from non-toxic kresion.
- Latex-free.

CODE	TYPE
WL431001	Without Vent



## Wellead Suction Tubing

- Suction tubing is made from non-toxic medical degree PVC.
- Without yankauer handle.
- Good transparency and tenacity.
- Latex-free.
- Batch marking for traceability.
- CE class IIA.
- Size: 1/4" x 360cm.

CODE	TYPE	SIZE
WL432002	Adult	1/4" x 360cm
WL432004	Paediatric	3/16" x 360cm


# OXYGEN MASKS & AIRWAYS



## Wellead Nasal Oxygen Cannula

- Made from clear, soft, non-toxic PVC.
- Latex-free.
- Batch marking for traceability.
- CE class IIA.

CODE	TYPE
WL351001	Adult



## Timesco Guedel Airways

- Latex-free.
- Individually wrapped.
- Rounded atraumatic edges.
- Bite blocks.
- Available from size 000 to size 4.

CODE	COLOUR	SIZE
TIME207	Pink	Size 000 - 40mm
TIME208	Blue	Size 00 - 50mm
TIME209	Black	Size 0 - 60mm
TIME210	White	Size 1 - 70mm
TIME211	Green	Size 2 - 80mm
TIME212	Yellow	Size 3 - 90mm
TIME213	Red	Size 4 - 100mm





PATIENT  
CARE

## PATIENT CARE

We offer a wide range of patient care products ensuring patients are provided with optimal care and comfort. The range includes bibs, cotton wool, sick bags dispensers and disposable slippers.



## Medicare Cotton Wool Roll

- 100% pure cotton.
- Maximum softness with minimal fibre shedding.
- Interleaved-easy to start rolls.
- Good structural characteristics.
- Pull string top for releasing after use.

CODE	SIZE
MD81	100g
MD83	500g



## BV Cotton Balls

The blend of premium quality cotton in BV Cotton Balls are kind to the skin.

- Non-sterile cotton balls.

CODE	COLOUR	PACK SIZE
BA444-001	White	100
BA444-002	White	250
BA444-003	White	500



## BV Non-Sterile Non-Woven Dry Wipes

- Made from 70% viscose 30% polyester.
- Absorbent.

CODE	SIZE	PACK SIZE
BA4561-300	30 x 35cm	100
BA4561-301	20 x 20cm	40



## Attends® Disposable Bib Regular 250's

Attends® disposable bibs are a perfect solution to keeping clothes dry and free from spillages, with the addition of an effective liquid absorbent top sheet to quickly absorb any spillages.

- The bib comes with adhesive fixation tape on the shoulders which is designed to keep the bib in place and ensure maximum security so you won't have the worry of the bib falling off during use.
- Polythene waterproof back sheet.
- Front pocket to reduce risk of spills.
- Absorbent top sheet.

CODE	SIZE	PACK SIZE
AT0117	67 x 35cm	250



## BV Impermeable Bib

- Blue polyethylene film laminated with white tissue paper.
- With fixation strips.
- Residue bag, completely waterproof.

CODE	SIZE	PACK SIZE
BA468-001/1	37 x 50cm	100



## BV Disposable Slippers

- Disposable slippers with closed front and anti-slip sole.
- Light-weight and breathable.
- One size fits all.

CODE	SIZE	COLOUR	PACK SIZE
BA468-005	27.5 x 10.5cm	White	800



### Medicare Sick Bags

Medicare disposable sick bags are easy to use, no-mess, no-fuss sick bags that lock away contaminants and odours for the comfort of the user and the safety of the people caring for them.

- Disposable.
- Compact.
- Easy and quick.
- No spill and odour-free.

CODE	PACK SIZE
MD594	50



### Sick Bags Dispenser

- Sick Bags Dispenser made of rigid white plastic which can be fixed to any kind of surface using screws.
- Capacity: 50 bags.

CODE
BA100-420



### Hygea Hands & Face Personal Washcloths

- For professional and hospital use for face and hands.
- Hypoallergenic solution which gently cleans and refreshes skin.
- Leaves skin clean and fresh.
- Significantly reduces bathing time and laundry cost.
- Pre-moistened.
- Alcohol and lanolin free to prevent irritation and cracking of skin.
- Easy to dispense.

CODE	SIZE	PACK SIZE
PD44	220 x 200mm	200





MATERNAL  
CARE

## MATERNAL CARE

Our maternal care range ensures the effective and comfortable care of maternity patients.





## Medicare Maternity Pads

Medicare Non-Sterile Maternity Pads are specially designed towels for new mothers. These soft and absorbent thick maternity pads cushion and protect, offering maximum comfort while featuring a unique slim-line design.

- Ultra absorbent.
- Slim design.
- No lateral leaks.
- Non-Sterile.

CODE	SIZE	PACK SIZE
MD4861	10 x 30cm	10



## Medicare Paediatric Urine Collector

- Capacity: 150ml.
- Display Box: 25 pcs.

CODE	CAPACITY	PACK SIZE
MD409	150ml	25



## BV Sterile Baby Blanket

- Sterile, soft, smooth heat retaining blanket.
- Made from 67% polyester and 33% viscose.
- Hypoallergenic and highly absorbent.

CODE	SIZE	PACK SIZE
BA454-001	90 x 80cm	25



## BV Non-Sterile Baby Blanket

- Non-Sterile, soft, smooth heat retaining blanket.
- Made from 67% polyester and 33% viscose.
- Hypoallergenic and highly absorbent.

CODE	SIZE	PACK SIZE
BA454-300	90 x 80cm	25



EMERGENCY

## EMERGENCY

Our range of emergency products provide essential items needed during an emergency such as stretchers, care bags and aspirator equipment.




## Evacuation Sheet EVS120

Evacuation sheet to be permanently in place under anti-pressure mattress 120cm wide. This fixes permanently to the bottom of a flexible foam mattress by means of strong elastics. It has safety belts to hold the patient on the mattress and pulling straps at each end. The combination of slippery sheet and mattress which will flex, will pass through narrow doorways and down stairs.

- The holding strap at the head is to hold back on stairways.
- No transferring of patient.

### CODE

HSP120



## Boardloc Spider Straps

This belts system allows a whole patient immobilisation. Thanks to its transversal sliding belts and velcro® straps, it can be used with any kind of patient and spinal board.

- Fully adjustable boardloc spider straps for spine board.
- Fits patients of varying sizes.
- Hook and loop enclosures provide secure connection to spine board.

### CODE

MEB01



## Head Immobiliser Block

- Head Immobiliser Block protects the patient's head and is suitable for use on most types of spineboards.
- Lightweight block made from closed cell foam.
- Easy to use and clean, impervious to blood and other body fluids.
- Large ear holes allow for monitoring of possible bleeding and/or fluid drainage.

### CODE

AD01/1



### Boscarol OB 3000 Suction Unit with Power Supply and Wall Bracket

- Medical Suction Unit OB 3000.
- Max. suction rate: 800 mbar (80 kPa).
- Max. suction flow: 33 litres/minute.
- Weight: < 2.8 kg.
- 1 litre liner version with disposable bag and FA version with autoclavable jar.
- Lithium polymer battery (solid-non organic).
- Battery status indicator with no need to turn on the device.
- Battery to be replaced by user.
- Works with 12 and 24 Vcc.
- Inner microprocessor for all data record and functioning tracing.
- Safety inspection reminder.
- Analogic fluorescent vacuum meter, with mmHg scale upon request.
- Compliant to all current norms, 93/42/CEE Directive and main reference norm ISO10079-1.

#### CODE

BOS38



### Boscarol OB 1000 Suction Unit with Power Supply and Wall Bracket

- Boscarol suction unit is powerful, compact and easy to use.
- Autonomy of approx 45 minutes continuous.
- Control-manometer from 800 mbar (80 kPa).
- Max. suction flow: 20 litres/minute.
- 1000ml jar with disposable liner.
- Direct operation through 12V C.C. of the vehicle via cable.
- AC adapter for car and recharging wall bracket.
- Weight: 3.6kg.
- Includes carry case.
- Conforms to ISO 10079-1-1999, EN1789 1999, classified on Class IIb and CE certified.
- 2 years Warranty.

#### CODE

BOS50





## Community Care Bag

- Its great capacity and simplicity allow you to transport everything required in home care.
- From basic to intuitive lines, it perfectly satisfies the professional's need in their visits.
- Waterproof main material for easy cleaning.
- Handle and shoulder strap for easy portability.
- Semi rigid structure of extreme lightness adaptable to the content.

CODE	SIZE	COLOUR
EB00018	38 x 28 x 18cm	Black



## TIMESCO Handheld Suction Unit

- Rescuer MVP Suction Pump Aspirator.
- Lightweight & portable.
- High vacuum.
- One hand operation.
- Simple maintenance.
- Multiple settings.
- Can be used neonates to large adults.

CODE
TIME110



## Boscarol Suction Liners

- Disposable 1000ml series suction liner bags.
- Works with Boscarol OB1000 Suction Unit (BOS50) and Boscarol OB2012 Suction Unit (BOS19).

CODE	CAPACITY	PACK SIZE
BOS50/2	1000ml	6

## TIMESCO BVM Resuscitator

The ion-air manual resuscitators are supplied complete with oxygen tubing, an oxygen reservoir bag and a suitably sized anatomical air cushion mask.

- Textured grip (non-slip).
- 360° swivel connector for greater versatility.
- Glued components (with exception of face mask).
- Adjustable air cushion face mask for maximum patient comfort.
- Transparent bag for complete internal visibility.
- Wrapped in a draw string bag.
- Single patient use.



CODE	SIZE	CAPACITY
TIME206	Infant Size 0	350ml
TIME1600	Adult Size 5	1600ml





INFECTION  
CONTROL

## INFECTION CONTROL YOU CAN TRUST

Our tried, tested and trusted gloves, hand sanitiser, bibs, face and disposable shoe covers, masks and wipes comply with European Quality Assurance Standards, ensuring safety and efficiency.



## Mediskin Latex Powdered Gloves

- Mediskin Lightly Powdered Latex Gloves are made from high quality natural rubber latex.
- Like other forms of barrier protection for the hands, lightly powdered latex gloves come in a range of sizes to accommodate different wearers.

CODE	SIZE	PACK SIZE
TG01	Small	100
TG02	Medium	100
TG03	Large	100
TG11	Ex. Large	100



## Mediskin Latex Powder-Free Gloves

- Mediskin Latex Powder-free Gloves are made from high quality natural rubber latex. Like other forms of barrier protection for the hands, powder-free latex gloves come in a range of sizes to accommodate different wearers.
- The latex examination gloves are 50% thicker than powdered gloves.

CODE	SIZE	PACK SIZE
TG04	Small	100
TG05	Medium	100
TG06	Large	100
TG12	Ex. Large	100



## Mediskin Nitrile Powder-Free Gloves

- Ideal for conducting medical examinations, diagnostic and therapeutic procedures to protect patient and user from cross contamination or infection.
- Non-Sterile Nitrile (Carboxylated Budadiene Acrylonitrile).

CODE	SIZE	PACK SIZE
TG60	Small	100
TG61	Medium	100
TG62	Large	100
TG71	Ex. Large	100



## Premier Glove Dispenser Rack

- Wall mountable glove dispensers.
- Manufactured from white plastic coated wire.
- Hard wearing and easy to clean.
- Supplied complete with fittings.

CODE
PM2135



## Medicare Hand Sanitiser

Probably the single most effective way of combating HCAI (Healthcare Associated Infections) is to improve hygiene in healthcare settings, in particular hand hygiene. Use anytime, any place without water or towel. Leaves hands feeling soft, smooth and refreshed without stickiness or residue.

- Anti-bacterial with 70% alcohol.
- Denturises the microorganisms cellular proteins.
- Kills 99.9% of most common germs.
- The gel contains moisturising properties.

CODE	SIZE
MD830	500ml



## Medipal Detergent Wipes

- 18 GSM material.
- 350 pieces per tub.
- High strength and ultra low linting Tritex material reduces risks of tearing or leaving surface fibre contamination.
- Low smear formula, providing effective cleaning power without leaving visible residue on glass and stainless steel.

CODE	SIZE	PACK SIZE
PAL535110	250 x 280mm	350



## Medipal Alcohol Disinfectant Wipes

- Kills 99% of germs, viruses, fungi and other microorganisms effectively and easily.
- The wipes are impregnated with 70% w/w isopropyl alcohol solution.
- These wipes are a convenient multi-tasking cleaning solution.

CODE	SIZE	PACK SIZE
PAL600110	190 x 210mm	200



## Medicare Disposable Face Mask

- 3 ply disposable non-woven face mask.
- Very low breathing resistance.
- Adaptable nose bar helps conform mask to your face.
- 95% bacterial filtration efficiency.
- Pleated folds allows user to expand the mask to cover the area from the nose to the chin.

CODE	SIZE	COLOUR	PACK SIZE
MD5221	175 x 95mm	Blue	50



## BV FFP2 Non-Sterile Face Masks

- Designed to protect against both solid and liquid aerosol.
- Reduces exposure to airborne particles.
- Produced with different non-woven plies, in order to satisfy its specific function. The different non-woven plies are ultrasonically welded.
- Filtration efficiency: >94%
- Individually wrapped.

CODE	PACK SIZE
BA465-102	20



## Disposable PE Overshoes

- Disposable PE overshoes.
- Multi-purpose shoe covers.
- The elastic keeps the cover securely in place.

CODE	PACK SIZE
W025	100



## Medicare Polyethene Aprons on a Roll

- Medicare Polythene Disposable Aprons are light, durable, disposable aprons with halter neck and waist ties.
- Easy non-stick removal of each apron.

CODE	SIZE	COLOUR	PACK SIZE
MD330/1	26 x 46"	White	200
MD331/1	26 x 46"	Red	200
MD333/1	26 x 46"	Green	200



## Premier Protect Barrier Gown Standard

This gown is manufactured from soft durable polyethylene material of 1 piece construction. Impervious to fluids, providing an extremely effective barrier against cross infection.

- Tear away halter neck allows easy removal of the gown.
- Independently tested as 100% effective against strike through by staphylococcus aureus (MRSA).
- Thumb loops prevent leakage between glove and sleeve.

CODE	COLOUR	PACK SIZE
PM5530	Blue	20
PM5550	Green	20



## Medicare Cleaning Cloths

Medicare Cleaning Cloths are made from 50% viscose and 50% polyester.

- Non-woven fabric.
- Colour: Blue.
- Size: 60 x 30cm.

CODE	PACK SIZE
MD700/1	50



## Disposable Single-Use Face Shield

- Transparent, shock resistant, economical and lightweight.
- Equipped with an elastic band and sponge headband for maximum comfort.
- Includes removable film to protect material from infection before use.
- FDA Class I medical device.
- Meets ANSI Z87.1-2015 eye protection standards for droplets and splash.

CODE	PACK SIZE
HPL01	100





## Beeswift Bergen Overspec Visitor Glasses

- Wraparound clear visitors spectacles.
- Polycarbonate lens.
- Ventilated side arms.
- Width: 157mm.
- Conforms to EN166 and ANSI Z87.
- Optical Class: 1.
- F - (Low impact energy).
- 2-1.2 - UV Filter.

### CODE

BS1035



## Delta Plus Milo Clear Safety Glasses

Featuring single lens, impact resistant polycarbonate lenses, the Delta Plus Milo Clear Safety Glasses use an ergonomic design to provide both excellent comfort and long-lasting protection.

They use a special integrated polycarbonate nose piece to ensure extra protection, while their flexible design looks to suit most users.

### CODE

SSPEC02



## HSE Metal Bin - White

- Fire Retardant Bin with anti-bacterial coating and slow, silent closing feature.
- Can be easily moved due to handle at rear.

CODE	SIZE	COLOUR
LM1000	20 litres	White
LM1001	50 litres	White
LM1002	70 litres	White



## HSE Metal Bin - Yellow

- Fire Retardant Bin with anti-bacterial coating and slow, silent closing feature.
- Can be easily moved due to handle at rear.

CODE	SIZE	COLOUR
LM1000Y	20 litres	Yellow
LM1001Y	50 litres	Yellow
LM1002Y	70 litres	Yellow



**MEDICINE  
STORAGE**

## **MEDICINE STORAGE**

Discover our range of pharmacy supplies and accessories such as tablet containers, ziplock bags and hygrometers. Our supplies are carefully designed to ensure the safe storage of medications, vaccines and critical supplies.



## Shoreline Solid Door Pharmacy Fridge

- Solid Door Medical Fridge suitable for the storage of pharmaceuticals and vaccines.
- RPSGB, NHS and RCVS compliant.
- Temperature range +2 to +8 degrees.
- Audio/visual alarm on temperature deviation.
- Automatic defrosting.
- Digital temperature display.
- Lockable door with 2 keys.

CODE	SIZE
SL45	45 litres
SL75	75 litres
SL161	141 litres
SL150	151 litres



## Shoreline Glass Door Pharmacy Fridge

- Solid Door Medical Fridge suitable for the storage of pharmaceuticals and vaccines.
- RPSGB, NHS, WHO, RCVS compliant.
- Temperature range +2 to +8 degrees.
- Forced-Air cooling for even temperature distribution.
- Digital temperature display.

CODE	SIZE
SL75G	75 litres
SL161G	141 litres

## Hygrometer with USB Data Logger and Software

Hygrometer with USB Data Logger and Software can measure either fridge temperature or room humidity and temperature.



- Memory for 32,000 readings.
- Dew point indication.
- Status indication.
- USB interface.
- User-selectable alarm.
- Analysis software.
- Multi-mode to start logging.
- Long battery life.
- Selectable measuring cycle: 2s, 5s, 10s, 30s, 1m, 5m, 10m, 30m, 1hr, 2hr, 3hr, 6hr, 12hr, 24hr.

CODE
ITS06



## Sample Containers

- Sterile Sample Container for both sputum and urine.
- Non-Spill cap.
- Plain label on each container.

CODE	SIZE	PACK SIZE
SP0001	60ml	60
SP0002/1	30ml with scoop	50
SP0003/1	30ml	50



## Zip Lock Bags

- Transparent resealable zip lock bag made from polyethylene.
- 2mm thickness.

CODE	SIZE	PACK SIZE
ZIP01	11 x 8.5cm	1000
ZIP02	13.5 x 10cm	1000
ZIP03	19 x 12.5cm	500
ZIP04	22.5 x 15cm	500



## Medicare Weekly AM/PM Push Button Pill Reminder

For travel, school, work, or daily life, nothing beats our Weekly Pill Organiser for keeping the right dose of your medications, vitamins, or supplements easily accessible.

- Extra large AM/PM compartments.
- Easy open pop-up cover.
- Easy grip button.

CODE

MD67585



## Medicare Deluxe Weekly Pill Organiser

Medicare Deluxe Weekly Pill Organiser dispenses four times a day: morning, noon, evening and bedtime.

- Durable, raised letters and braille markings for the sight-impaired.
- Removable daily sections.
- Pull out medication scheduling card slides easily inside planner.
- Size: 25 x 17 x 3.5cm.

CODE

MD400407



## Medicare Weekly Pill Organiser

Medicare Weekly Pill Organiser dispenses four times a day: morning, noon, evening and bedtime.

- Seven removable daily compartments for travelling.
- Durable.
- Lids click securely closed.
- Size: 17.5 x 15 x 2.5cm.

CODE

MD128BW



## Medicare Weekly AM/PM Contoured Pill Reminder

Medicare Weekly AM/PM Contoured Pill Reminder with separate AM and PM compartments.

- Contoured bottom rounded for easy pill removal.

CODE

MD67585





## Medicare Round Weekly Pill Reminder

High quality, compact 7 day medication organiser.

- Durable, raised letters and Braille markings for the sight-impaired.
- Lids click securely closed.
- Unique, palm-sized shape.

### CODE

MD104C



## Medicare Weekly Push Button Pill Reminder

7 Day Push Button Pill Reminder.

- One Compartment per day.
- Easy open, pop-up cover.

### CODE

MD67577



## Medicare Small 7 Day Pill Reminder

Durable, raised letters and braille markings for the sight-impaired

- Sturdy hinges and easy-open lids with lips for gripping.
- Convenient size for out of home use.
- Curved bottom for easy pill removal.
- Lids click securely closed.

### CODE

MD100B



## Medicare Large 7 Day Pill Reminder

A sturdy design, with high visibility with each compartment clearly labelled for each day of the week.

- Braille markings for the sight-impaired.
- Convenient size for out of home use.
- Curved bottom for easy pill removal.

### CODE

MD300B



## Medicare One Day Pill Organiser

The Medicare One Day Pill Organiser dispenses four times a day: morning, noon, evening and bedtime.

- Helps with pill compliance.
- Easy to open and clicks safely back in to place.

### CODE

MD188B



## Medicare Daily AM/PM Push Button Pill Planner

The AM-PM Push Button Pill Planner offers easy access to daily medication.

- Easy access to compartments without having to use two hands.
- Great for those with dexterity issues.
- Compact size fits easily in pocket or purse.

### CODE

MD67579



## Medicare Daily AM/PM Pill Planner

Daily AM/PM Pill reminder.

- Handy Pocket size.
- Rounded scoop bottom for easy pill removal.
- Each top clearly marked AM and PM for the visually impaired.

### CODE

MD67433



## Medicare Pill Splitter

Consists of a stainless steel blade which splits all size tablets.

- Built-in compartments.
- Snap-shut door.
- Small and easy to carry.
- Dishwasher safe.

### CODE

MDPS12E



## Medicare Pill Splitter & Crusher

The combo splitter and crusher is a handy all-in-one tool.

- Stainless steel blade in lid for splitting pills.
- Built-in pill storage area.
- Bottom serves as a drinking cup.
- Separates easily and crushes pills with one twist.

### CODE

MD400087














CONTINENCE  
CARE

## CONTINENCE CARE

Browse our complete range of high-performance product solutions for continence care. Fleming Medical are a distributor of advanced Attends® and MoliCare® products.

## Attends® Absorbency Guide

Light		Moderate	Heavy		
0		2			
1		3			
		4		6	
		5		7	
				8	
				9	
				10	



### Attends® Pull-Ons

Attends® Pull-On are absorbent disposable pants designed to manage heavy bladder and/or bowel incontinence (urinary/faecal). These pants look like regular underwear to help you feel confident throughout the day, and can be worn without any other underwear required. The pants have a flexible shape which fit close to the body and an absorbent core to help reduce the risk of leakage.

- Can easily be pulled up and down like normal underwear.
- Superior fit and feel.
- 100% Breathable for All-around Air Comfort.
- Maxi Comfort Technology. Full panel elastics with ultrasonic bonding.
- The powerful absorbent core locks in any liquid and odour.
- Soft and skin-friendly: The skin tolerability of Attends® Pull-Ons is dermatologically approved by proDERM Institute for Applied Dermatological Research.

CODE	LEVEL	SIZE	WAIST/HIP SIZE	ABSORBENCY	PACK SIZE
AT211748	Level 4 	Small	56-75cm	850ml	22
AT211786	Level 4 	Medium	75-110cm	1049ml	22
AT211823	Level 4 	Large	100-140cm	1095ml	22
AT211953	Level 6 	XX Small	38-55 cm	1400ml	21
AT211977	Level 6 	X Small	45-60cm	1400ml	22
AT211762	Level 6 	Small	56-75cm	1519ml	18
AT211809	Level 6 	Medium	75-110cm	1620ml	18
AT211847	Level 6 	Large	100-160cm	2264ml	18
AT211885	Level 6 	X Large	120-175cm	2069ml	16
AT0028	Level 8 	Medium	80-110cm	2529ml	16
AT0024	Level 8 	Large	100-160cm	2529ml	16





### Attends® Cover-Dri Underpad Super

Attends® Cover Dri underpads are designed to provide added surface protection for beds, mattresses, chairs and wheelchairs. The product has a unique hexagonal pattern with super absorbent materials that quickly lock urine into the pad. Attends® Cover Dri has a soft, skin-friendly top sheet that provides extra dryness and comfort.

- Has a unique, hexagonal design that offers improved dryness and comfort for incontinence sufferers.
- Cover-Dri underpads can be used for short, medium and long term use.
- Channels between the hexagonal cushions do not intercede, reducing the risk of 'pooling'.

CODE	SIZE	PACK SIZE
AT0058	90 x 60cm	50



### Attends® Cover-Dri Underpad Plus

Attends® Cover Dri underpads are designed to provide added surface protection for beds, mattresses, chairs and wheelchairs. The product has a unique hexagonal pattern with super absorbent materials that quickly lock urine into the pad. Attends® Cover Dri has a soft, skin-friendly top sheet that provides extra dryness and comfort.

- Has a unique, hexagonal design that offers improved dryness and comfort for incontinence sufferers.
- Cover-dri plus underpad is designed for medium term use.

CODE	SIZE	PACK SIZE
AT0050	60 x 60cm	50



### Attends® Cover-Dri Plus











Attends® Cover Dri underpads are designed to provide added surface protection for beds, mattresses, chairs and wheelchairs. The product has a unique hexagonal pattern with super absorbent materials that quickly lock urine into the pad. Attends® Cover Dri has a soft, skin-friendly top sheet that provides extra dryness and comfort.

- Hexagon design delivers improved protection and faster absorbency.
- Strong, anti-slip backsheet, ensuring the pad will stay in place.

CODE	SIZE	PACK SIZE
AT0078	90 x 60cm	50



## MoliCare® Drop Level Guide


Slight	Moderate	Severe	Very Severe
1 	4 	7 	9 
2 	5 	8 	10 
3 	6 		



### MoliCare® Premium Mobile Level 5 & 6

MoliCare® Premium Mobile 5 & 6 Drop Level pull-up pants look and feel like regular underwear through anatomical shape, soft material and elastics for a comfortable fit.

- The absorbent core neutralises odour and holds liquids quickly for a dry feeling.
- Easy tear-open side seams mean they can be easily pulled on and off without the need to open and close any tapes or additional help.



CODE	LEVEL	SIZE	SIZE	ABSORBENCY	PACK SIZE
MOL01	Level 5 	Large	60 x 90cm	680ml	14
MOL02	Level 6 	Medium	80 x 120cm	680ml	14



### MoliCare® Premium Bed Mat

MoliCare® Bed Mat products provide extra protection for various surfaces and the Rectangular pads are the classic protection alternative for pads.

- Hygienic Care: Single-use protection.
- Absorbs liquid and provides extra protection for various surfaces.
- Dermatologically tested and skin-friendly.
- Grip edging allows for easier unfolding of bed mat.

CODE	LEVEL	SIZE	PACK SIZE
MOL05	Level 5 	60 x 90cm	100
MOL03	Level 7 	60 x 90cm	30

## MEDICAL CONSUMABLES

We offer the highest quality medical consumables for patient care, including airways, oxygen masks, nasal cannula, infusion and sterilisation products, procedure packs, forceps and clamps, sharps, blades and more.



## Disposable Pen Torches

- Disposable pen torch is ideal for general practice and everyday medical applications.
- Pupil gauge printed on torches to aid diagnosis.
- Convenient pen clip also switches torch on when pushed against barrel.

CODE	COLOUR	PACK SIZE
ML37	White	6



## BV Gallipots 60ml

- Gallipot containers are made from 100% polypropylene.
- Sterilisable by steam or gas (ETO).
- Graduated in 10ml increments.

CODE	CAPACITY	PACK SIZE
BA489-304/3	60ml	1,200



## Premier Tongue Depressors (Wooden)

- Non-sterile, flat, wooden spatula.
- Smooth, safety rounded ends.
- Produced and packed under clean conditions.
- Each box is cellophane wrapped to protect from dust.

CODE	PACK SIZE
PM1604	100



## BV Plastic Adult Curved Tongue Depressors

- Plastic Adult Tongue Depressors are used for oral patient examination and are also useful for stirring and mixing medication and spreading ointments.
- Made from medical grade polystyrene.
- Clear plastic smooth edges.

CODE	SIZE	PACK SIZE
BA469-103	15 x 2cm	100



## Medicare Foam Cleaning Sticks

- Single use.
- Ideal for use in oral care.
- Used to lubricate the mouth and lips.

CODE	SIZE	PACK SIZE
MD410	150 x 4cm	5



## BV Crepe Paper

- Medical grade paper - 60g/sq.m.
- For sterilisation by steam, ETO - ethylene oxide (gas) and irradiation.

CODE	SIZE	COLOUR	PACK SIZE
BA450-006/1	50 x 50cm	Blue	500



## BV Drawsheets Non-Sterile

- Polyethylene blue film with white paper tissue.
- Resistant, with an impermeable side and an absorbent side.
- Transport carton.

CODE	SIZE	PACK SIZE
BA455-102	150 x 80cm	420



## BV Drawsheet

- Laminated polyethylene film with tissue paper.
- Resistant, with an impermeable side and an absorbent side.

CODE	SIZE	PACK SIZE
BA455-303	150 x 80cm	140
BA455-304	180 x 80cm	120





## BV Examination Couch Covers

- Non-woven polypropylene.
- Loose blanket.

CODE	SIZE	COLOUR	PACK SIZE
BA4541-003	75 x 210cm	White	48



## BV Paper Exam Table Covers

- Pre-cut every 38cm.
- 2 ply.

CODE	SIZE	COLOUR	PACK SIZE
BA4541-101	48cm x 50m	White	130 sheets



## BV Mattress Covers

- Polyethylene film.
- Thickness: 30 microns.
- Impermeable to fluids and with elastics to adjust to the mattress.

CODE	SIZE	COLOUR	PACK SIZE
BA4541-104	210 x 90 x 20cm	Blue	10



## Heavy-Weight Body Bag

- Adult Size.
- Centre Zip.
- PVC body bag.

CODE	SIZE	COLOUR	PACK SIZE
BB01	228 x 75cm	White	20



## NOTES

## NOTES





## Here to Help

For support and further information

**FLEMING MEDICAL LTD.**

Corcanree Business Park, Dock Road,  
Limerick, V94 K6W4, Ireland

T: 1800 307777

E: [healthcare@flemingmedical.ie](mailto:healthcare@flemingmedical.ie)



EN ISO 13485:2016

**FLEMINGMEDICAL.IE**



[fleming-medical-ltd](https://www.linkedin.com/company/fleming-medical-ltd)



[@FlemingMedical](https://www.facebook.com/FlemingMedical)



[@FlemingMedical](https://twitter.com/FlemingMedical)