



Medicare Burn Gel with Aloe Vera

Hydrogel with soothing Aloe Vera, promoting faster wound healing



FACT BOX

Medicare Burn Gel MD1220

Product Category: Wound and skin care, First Aid, OTC.

Composition: Purified water, Carbomer, Preservatives, Aloe Vera, Sodium Hydroxide.

Product Approvals: Registered according to MDD (Medical Device Directive, 93/42/EEC, Class IIa)

Shelf Life: 3 years from production date

Packaging: Standard can size is 100ml.

Precautions: Do not use on deeper second degree burns or other types of burns such as chemical burns, phototoxic reactions, drug induced skin reactions or skin damages induced by electrical shock.

Medicare Burn Gel with Aloe Vera is a water based gel for the first aid treatment of sunburns, minor burns, and scalds.

The gel works by forming a transparent protective film as protective barrier that cools and soothes the skin to help reduce pain and discomfort.

It is effortless to apply while also being water soluble and is easily washed off.

Bag-on-valve technology allows you to use every last drop with a longer shelf life. No pumping motion needed.

FEATURES

- Cooling effect on burns and pain relief.
- Transparent protective gel.
- Non-contact, easy and comfortable application.
- Can be used with the canister in any position.

DIRECTIONS:

Apply the gel abundantly on the affected area immediately following a thermal injury. Apply from a close distance.

RECOMMENDED USE:

Can be used by children and adults. Store between 15°C (59°F) and 30°C (86°F).